# SOLUTIONS FOR CHROMATOGRAPHY AND PURIFICATION









Founded in 1995, SiliCycle is specialized in the development, manufacturing and commercialization of high value silica gels and specialty products for chromatography, purification and synthesis.

# **Solutions for Purification and Chromatography**

Silia <i>Flash</i> ® and Silia <i>Sphere</i> ™ – Bulk Silicas for Chromatography	3
Silia <mark>Bond</mark> ® – Chromatographic Phases	13
SiliaSep™ – Flash Cartridges	19
Silia <i>Plate</i> ™ – Thin Layer Chromatography Plates	31
Resource Center	47





SiliCycle is your partner of choice for your purification and chromatography needs! Recognized as one of the leaders with an excellent quality silica gel, SiliCycle offers a wide range of products available in two different shapes:

- Silia*Flash*® Irregular silicas
- SiliaSphere<sup>™</sup> PC Spherical silicas



# SiliCycle: Silica Expert

With pore diameters ranging from 30 to 1,000 Ångström (Å) and particle sizes up to 1,200 microns ( $\mu m$ ), SiliCycle offers products to meet all your requirements. We offer one of the most reliable portfolios for flash and gravity grades for low to medium-high pressure. Our silica gels are ideal for preparative chromatography, from laboratory to pilot-plant processes and production scale.

Features and Benefits of SiliaFlash & SiliaSphere PC					
Features	Benefits				
High purity silica gels	No contamination, consistency, reliability, reproducibility				
Low level of fines	No contamination, lower back-pressure, good separation				
Narrow particle and pore size distributions	Optimal separation and resolution				
Batch-to-batch, year-to-year consistency	Reliable chromatography				
Neutral pH	Wide range of products can be purified, even acid sensitive ones				
Low metal content and controlled water content	Symmetrical peaks without tailing				
High mechanical stability	Can be used under high pressures without surface abrasion				
High surface area and density	Greater loading capacity, enabling more silica for the same volume Solvent economy (smaller dead volume)				
Availability in bulk quantities	In stock for fast delivery				

Together, all these benefits mean optimal and reproducible separation power, saving you time and money.

# SiliCycle, the Silica Supplier for Every Need

Each year, SiliCycle manufactures hundreds of tons of silica for a broad range of chromatography applications. All our products are manufactured under tightly controlled manufacturing processes and a stringent quality control ensures the highest quality.

Be confident in scaling-up your processes with our silica gels.

With SiliCycle, scale-up from laboratory to production without limitations!







Enjoy a virtual tour of SiliCycle's facility





Path straight

through the column

# Two Shapes Available: Irregular and Spherical

The quality of a silica gel is extremely important when you are using it for chromatography purposes, particularly when dealing with difficult separations of valuable compounds. You need to be confident about your recoveries.

In chromatography, there are at least three physical properties that will influence your separation and that you need to consider when choosing your silica gel:

- Particle shape (irregular or spherical)
- Particle size distribution (tight or large)
- · Pore diameter (surface area)

These caracteristics will directly influence crucial parameters involved in a successful chromatography:

- Resolution (efficiency of separation and final purity)
- Retention (which allows separation)
- Capacity (maximal sample quantity and final recovery / yield)
- Back-pressure (speed and pumps related issues)

At SiliCycle, we ensure consistency, reliability and reproducibility.

Our expertise and strong knowledge has been developed over many years of helping our customers find the best solutions to their particular needs.

# How to Choose Between Silia Flash Irregular and Silia Sphere PC Spherical Gels?

Irregular silica gels are traditional in flash or gravity chromatography and have always been a spontaneous choice for preparative chromatography. Nowadays, spherical particles are used increasingly.

Cost is very important in preparative and process chromatography, and the use of monodisperse spherical particles with narrow particle size distribution is more expensive. It is possible in this case to use irregular silica but the separation may not provide the desired results. For these situations, SiliCycle has developed a more affordable class of spherical particles for preparative chromatography: SiliaSphere PC.

Advantages of using Silia Sphere PC materials over standard irregular silica gels include the following:

- · Increased efficiency of the eluent's flow characteristics
- Higher resolution
- Ease of packing / better packing reproducibility
- · Higher mechanical stability

# SiliaSphere PC: Truly Spherical

Silica gel quality varies greatly between manufacturers. Even when advertised as being "spherical" this may not be the case. Please discover on next page a quick comparison of Scanning Electron Microscopy (SEM) pictures between SiliCycle Silia. Sphere PC and the competition.



# Silia Flash and Silia Sphere PC Characteristics


The importance of the particle and pore size distributions varies depending on the type of chromatography being done

Importance of Tight Distributions in Chromatography					
Tight Particle Size Distribution	Tight Pore Size Distribution				
Greater column performance and separation	Surface area (Presence of bigger pore size leads to lower surface availability)				
Tighter peaks and better peak shape	Optimal peak shape (Presence of smaller pore size leads to peak tailing)				
Better column packing, easier to pack					
No preferential pathways (channeling)	No molecule sequestration due to fluid diffusion inside pores				
Faster flow rate with lower back-pressure	No molecule sequestration due to iluid diffusion inside pores				
Time and solvent savings					

### Scanning Electron Microscopy (SEM) comparison of two IRREGULAR silica gels 40 - 63 µm, 60 Å







SiliCycle

Competitor

SiliCycle

### **Effects of Homogeneous vs Uneven Packing**

The connection between particle size distribution and column performance is very simple. When the distribution is broad, the packing is uneven. Some parts are composed of only large particles where the solvent will flow fast and meet little resistance, and there are sections composed of small particles where the solvent flows slowly and meets great resistance.

As a result, the solvent will take the path of least resistance through the column and flow around the pockets of small particles instead of straight through the column.

This uneven flow greatly affects the separation because the compounds will have different retention times depending on their flow path. As they exit the column, the compounds will give broad and poorly separated peaks.


The figure on the right illustrates the effect of a wide particle size distribution versus a narrow one. Narrower distribution gives a more homogenous packing and thus more concentrated fractions. And, by reducing solvent consumption, the process will be more cost-efficient.




Narrow Particle Size Distribution

Scanning Electron Microscopy (SEM) comparison


of two SPHERICAL silica gels 50 µm, 60 Å






Competitive Gel


SiliaFlash Silia Sphere PC





Preferential


path chosen





# **High Purity Silica Gels**

You can be sure of the outstanding quality of SiliCycle's silica gels because of the closely controlled manufacturing conditions. Our tight control of every manufacturing process step allows reproducible results (chemical, physical and structural) as well as ensuring the same chromatographic selectivity. Hence, Silia Flash and Silia Sphere PC are suitable for validated chromatographic processes.





#### **Stable Water Level Content**

Water level of silica gel affects the selectivity of the silica. Silia Flash and Silia Sphere PC have generally a water content between 2 to 6 %. This is advantageous for you since other products have a water variation from 2 to 15 % depending on the manufacturer. SiliCycle can also adjust the water level upon request.

#### Neutral pH

Our silicas are pH-adjusted between 6 and 8 to be safely used in the separation of a wide range of products (a neutral pH is needed to separate pH-sensitive compounds). Once again, this is advantageous when compared to many gels on the market that are much more

#### **Low Trace Metal Content**

Silica, depending on its method of manufacturing, contains a certain amount of various metals. This can, in turn, affect the quality of the separation. Aluminum, iron and lead are particularly problematic because they cause peak tailing. SiliCycle's proprietary technology generates a silica gel with the lowest trace metal content on the market. This ensures you will get optimal performance from your chromatography. Tight control of metals in every batch also improves your reproducibility and reduces risks of interaction between metals and desired compounds.

Typical Metal	Content Con	nparison for 40 - 63	μm, 60 Å Silica Gel	ls (mg/kg)
Metals		SiliCycle F60 R10030B	Manufacturer A	Manufacturer B
Aluminum	Al	33	262	280
Barium	Ва	9	60	33
Calcium	Ca	336	1,150	502
Chromium	mium Cr		0.6	0.4
Copper	Cu	0.2	0.2	0.2
Iron	Fe	32	75	41
Lead	Pb	0.41	5.3	0.95
Magnesium	Mg	61	149	104
Nickel	Ni	0.4	0.5	0.5
Silver	Ag	0.09	0.29	0.19
Sodium	Na	466	945	585
Tin	Sn	0.2	0.2	0.1
Titanium	Ti	147	250	179
Zirconium	Zr	32	75	56

# Silia*Flash* Irregular Silica Gels

### One of the Tightest Particle Size Distribution on the Market

#### **Particle Size Distributions' Disparity**

When selecting a silica gel, chemists need to take into account that not all 40 - 63 µm gels are the same.


In this example, the figure on the right shows the distribution curves of SiliCycle's SiliaFlash gel (PN: R10030B) compared to other manufacturers of flash silica gels of same particle sizes. All products were sold as 40 - 63 µm 60 Å gels.

As you can observe, SiliCycle's gel has a mean of 90 % of the particles in the nominal range compared to maximum 80 % for the competitor gels. The higher the curve, the tighter the particle size distribution.

### Importance of the Absence of Fines

In chromatography, fine particles (small particles under 10 microns) increase back-pressure and can result in clogging. which is particularly dangerous when using glass columns. Fines can also pass through filters and contaminate final products. The lack of fines gives a more regular, stable and reproducible chromatography bed and a faster and more even flow rate for better separation.

The zoomed part of the figure shows that our most popular silica gel, Silia Flash 40 - 63 microns 60 Å, has total absence of fines unlike the six competitor gels analyzed.




# **Particle Size Analysis Methods**


### Laser Diffraction (Malvern Analysis)

Typically used for particle sizes below 40 microns. Particle size distribution is reported in term of D10, D50 (average, mean) and D90. Some manufacturers also mention the ratio of D90/D10.

#### Sieving

Usually for particle sizes over 40 microns. Particle size distribution is reported in percentage of undersized and oversized.









### **Two Different Grades for Different Needs**

Over the years, SiliCycle has developed two different grades ("Superior" and "Standard") for the two most popular irregular gels used in the industry:


- 40 63 μm, 60 Å
- 60 200 µm, 60 Å

Those two grades of each gel are available to address all our customers' requirements, depending on their applications, areas of research, budgets and so on.

### 40 - 63 μm, 60 Å Gels: Silia Flash F60 (R10030B) VS Silia Flash P60 (R12030B)

Both compare favorably with the overall industry average of a  $40 - 63 \mu m$  distribution, and each grade offers its own particle size distribution profile.

Tw	ım, 60 Å Gel	
Grade	Superior Grade	Standard Grade
Name	F60	P60
PN	R10030B	R12030B
Particle Size	40 - 63 μm	40 - 63 μm
Pore Diameter	60 Å	60 Å
Particularities	Extra step to reduce metal content to minimum level     Tighter particle size distribution     Fines have been removed	Fines have been removed     Lower price




The figure on the right shows F60's tighter particle size distribution and the absence of fines for both gels.

### 60 - 200 μm, 60 Å Gels: Silia*Flash* G60 (*R10040B*) VS Silia*Flash* GE60 (*R10140B*)

Each grade offers its own particle size distribution profile.

Two Different Grades of 60 - 200 µm, 60 Å Gel								
Grade	Superior Grade	Standard Grade						
Name	G60	GE60						
PN	R10040B	R10140B						
Particle Size	60 - 200 μm	60 - 200 μm						
Pore Diameter	60 Å	60 Å						
Particularities	Extra step to reduce metal content to minimum level     Tighter particle size distribution     Fines have been reduced to minimal level	Fines have been reduced to minimal level     Lower price						



The figure on the right shows G60's tighter particle size distribution.

# Typical Metal Content Comparison Between SiliCycle's Five Most Popular Gels

$\mathcal{L}$		Typical Meta	l Content of Most Pop	oular Irregular Silicas			
Product		F60	P60	Acid Washed	G60	GE60	
Product Number		umber R10030B R12030B R10530B			R10040B	R10140B	
Particle Size			40 - 63 μm	60 - 2	00 μm		
Pore Diameter			60 Å		60	Å	
Metal ( <i>mg/kg</i> )							
Aluminum Al		Al < 200 < 1,000 < 70		< 70	< 350	< 900	
Antimony	Sb		< 0.2		< (	0.2	
Arsenic	Ar		< 1		<	1	
Barium	Ва	< 40	< 40	< 5	<	40	
Beryllium	Ве		< 0.1		< (	0.1	
Bismuth	Bi		< 1		<	1	
Cadmium	Cd		< 0.01		< 0	.01	
Calcium	Ca	< 200	< 500	< 10	< 250	< 500	
Chromium	Cr		< 1	<	1		
Cobalt	Co		< 0.1		< (	0.1	
Copper	Cu		< 1		<	<1	
Iron	Fe	< 75	< 350	< 10	< 75	< 350	
Lead	Pb		< 1		<	1	
Lithium	Li		< 0.1		< (	0.1	
Magnesium	Mg	< 150	< 250	< 10	< 100	< 150	
Manganese	Mn	< 1	< 2	< 1	<	1	
Molybdenum	Мо		< 0.1		< (	< 0.1	
Nickel	Ni		< 1		<	1	
Potassium	К	< 500	< 30	< 2	< 750	< 30	
Rubidium	Rb		< 0.2		< (	0.2	
Selenium	Se		< 1		<	1	
Silver	Ag		< 0.1		< (	0.1	
Sodium	Na	< 150	< 1,500	< 15	< 150	< 1,500	
Strontium	Sr	< 4	< 15	<1	< 4	< 15	
Tellurium	Те		< 0.1		< (	0.1	
Thallium	TI		< 0.1		< (	0.1	
Tin	Sn	< 0.4	< 0.4	< 0.2	< 0.4		
Titanium	Ti	< 200	< 250	< 90	< 2	250	
Uranium	U		< 0.1		< (	0.1	
Vanadium	V		< 1		<	1	
Zinc	Zn		< 1		<	1	






Solutions for Chromatography and Purification

# **Irregular Silica Selection Guide**

Selecting the most appropriate sorbent for any given application can be difficult. To help you choose the right pore diameter and particle size, simply follow the two pathways to select the most suitable sorbent.

### **Selecting Pore Diameter**




### **Selecting Particle Size**

Generated Pressure (GP) or Selectivity (α)



Selectivity ( $\alpha$ ) is measured by the retention factor ratio between two similar compounds:  $\alpha$  =



# A Particle Size for Each Application

Y.		Most Popular Particle Size Applications								
Particle Size Irregular Particles	Distribution Spherical Particles	<b>Applications</b>								
Particles for Preparative	TLC Plates									
From 0 to 20 µm	-	<ul> <li>Contains neither binder (<i>organic or inorganic</i>) nor UV indicator (F<sub>254</sub>)</li> <li>Can also be used in flash chromatography if higher resolution is required (<i>higher back-pressure</i>)</li> </ul>								
Particles for Difficult Sep	arations									
From 10 to 45 μm From 15 to 45 μm • High-resolution silica for difficult separations (similar polarities)										
Particles for Flash Chron	natography									
40 - 63 μm	From 40 to 75 μm	<ul> <li>Chromatography types:</li></ul>								
60 - 120 μm	From 60 to 150 μm	Alternative to 40 - 63 μm silica for faster flow rate with lower pressure								
Particles for Column (or	Gravity) Chromatography									
From 60 to 200 μm	From 75 to 250 μm	Most economical silica for open column chromatography (gravity)     Suitable for very dirty purification     Easier to handle								
From 120 to 200 μm	From 100 to 200 μm	Silica for standard open column chromatography     Narrow particle size distribution enables uniform packing     Suitable for mass overload purification								
Other Application										
From 200 to 1,000 µm	From 200 to 500 μm	Silica for plugs								







Solutions for Chromatography and Purification

# Silia Flash & Silia Sphere PC Ordering Information

This is only an overview of gels we can provide. Please contact us if you are looking for a different product: <a href="mailto:support@silicycle.com">support@silicycle.com</a>. Available formats: from 1 kg to 25 kg, even up to multi-ton scale!

SiliaFla	ısh Irregular Si	lica Gels Por	tfolio	$\mathcal{C}$
Product Number	Partic µm	le Size mesh	Pore Diameter (Å)	Proc
R10137L	75 - 150	100 - 200	30	S100
R10130A	40 - 63	230 - 400		S100
R10150A	60 - 120	325 - 625		S100
R10140A	60 - 200	70 - 230	40	S100
R10160A	120 - 200	70 - 125	40	S100
R10170A	200 - 500	35 - 70		S100
R10180A	500 - 1,000	18 - 35		S100
R10117B	15 - 40	*		S100
R10023B	20 - 45	*		S100
R10030B (F60)			]	S100
R12030B ( <i>P60</i> )	40 - 63	230 - 400		S100
R10530B (Acid-Washed	n			S100
R10050B	60 - 120 325 - 625	1	S100	
R10040B (G60)	00 000	70 000	60	S100
R10140B (GE60)	60 - 200	70 - 230		S100
R10137B	75 - 150	100 - 200	1	S100
R10160B	120 - 200	70 - 125	1	S100
R10165B	150 - 250	60 - 100	1	S100
R10170B	200 - 500	35 - 70	1	S100
R10180B	500 - 1,000	18 - 35	1	S100
R10130D	40 - 63	230 - 400		S100
R10140D	60 - 200	70 - 230	1	S100
R10170D	200 - 500	35 - 70	90	S100
R10180D	500 - 1,000	18 - 35	1	S100
R10181D	800 - 1,200	16 - 22	1	S100
R10130H	40 - 63	230 - 400		S100
R10150H	60 - 120	325 - 625	1	S100
R10140H	75 - 250	60 - 200	1	S100
R10160H	120 - 200	70 - 125	150	S100
R10170H	200 - 500	35 - 70	1	* Mesh
R10180H	500 - 1,000	<del></del>		
R10181H	800 - 1,200	16 - 22	1	
R10130M	40 - 63	230 - 400		
D40440N4	60, 000	70 000	1 200	

70 - 230

35 - 70

60 - 200

200 - 500

SiliaSphe	ere PC Spherica	al Silica Gels	Portfolio		
Dundret Namelan	Partic	le Size	Dava Diameter (8)		
Product Number	μm	mesh	Pore Diameter (Å		
S10095W-A	25	*	50		
S10030B-A	50	300			
S10040B-A	300	150	- 60		
S10020C	20 - 45	*			
S10030C	40 - 75	200 - 400	70		
S10040C	75 - 200	70 - 200	- 70		
S10070C	200 - 500	35 - 70			
S10095D-A	25	*	90		
S10020E	20 - 45	*			
S10030E	40 - 75	200 - 400	100		
S10040E	75 - 200	70 - 200	100		
S10070E	200 - 500	35 - 70			
S10027G-A	50	300	120		
S10020M	20 - 45	*			
S10030M	40 - 75	200 - 400			
S10040M	75 - 200	70 - 200	300		
S10070M	200 - 500	35 - 70			
S10020P	20 - 45	*			
S10030P	40 - 75	200 - 400			
S10040P	75 - 200	70 - 200	500		
S10070P	200 - 500	35 - 70			
S10020S	20 - 45	*			
S10030S	40 - 75	200 - 400	000		
S10040S	75 - 200	70 - 200	- 800		
S10070S	200 - 500	35 - 70	1		
S10020T	20 - 45	*			
S10030T	40 - 75	200 - 400	1 000		
S10040T	75 - 200	70 - 200	1,000		
S10070T	200 - 500	35 - 70	1		

equivalent too small to exist as real screen size

### R10530B: Acid-washed Silia*Flash* 40 - 63 µm, 60 Å irregular silica gel for extra purity

300

This product gel has been developed to ensure a pH-controlled media with even lower levels of trace metal contaminants and maximal purity.

# SILICYCLE (4)

R10140M

R10170M

# **Chromatographic Phases**

Thanks to its high mechanical resistance, silica is the most widely used media in chromatography. With SiliaBond irregular silica gels, SiliCycle offers a large range of solutions for low pressure chromatography, to help cover many kinds of purification.

We quarantee quality and stability of our phases: no fines will appear when packing the media. Our gels will give you excellent performance and lifetime!

# Silia Bond Phases for Low Pressure Chromatography

For all our listed SiliaBond sorbents, particle size is 40 - 63 µm and pore diameter is 60 Å. Contact us if you need a different particle size or pore diameter: support@silicycle.com.

All functionalized Silia*Bond* sorbents are available in bulk but also pre-packed in Silia*Sep* flash cartridges and Silia*Prep* SPE cartridges.

### **Reversed-Phases**

In reversed-phase chromatography, the packing material is always hydrophobic (non polar) while the mobile phase is polar. The more hydrophobic the packing material, the more retention of non polar analytes.

Usual reversed-phases are standard alkyl chains grafted on silica (C18, C8, C4, C1) and cyclic or aromatic functions (Phenyl, Pentafluorophenyl).

Important parameters to keep in mind in reversed-phase chromatography:

- Carbon load (% C) will give the relative hydrophobicity of the packing media. Most of the time, it varies between 5 % and 17 %.
- Endcapping: when functionalizing silica, it is impossible to react with all available silanol groups (free -OH groups on the silica surface). But these free silanols are acidic and will react with basic compounds, so we endcap them with a capping agent to avoid non-specific bindings.


### **Normal Phases**

In normal phase chromatography, the packing material is always polar while the mobile phase is non polar. The interactions between analytes and sorbent mainly take place on the highly polar silanols of the silica gel surface. Some hydrogen bonds can also happen on polar functionalized groups.

Usual normal phases are ungrafted silica, polar functions (amine, cyano and diol) or alternative adsorbents (Alumina and Florisil for example).

# Ion Exchange Phases


In ion exchange chromatography, both silica support and analytes must be ionized. If the stationary phase (packing material) is positively charged, anionic analytes only will retain (these phases are called WAX and SAX). And in the contrary if the stationary phase is negatively charged, cationic analytes only will retain (these phases are called WCX and SCX). Hence, pH of the mobile phase is of crucial importance and needs to be chosen carefully, so that both functions are charged:





# **Sorbent Selection Chart**

SiliCycle offers a wide range of SiliaBond sorbents to cover many kinds of purification. The following chart will guide you for the selection of the appropriated sorbent, based on the characteristics of the sample to purify.



# **Reversed & Normal Phases Portfolio**

Available formats: from 5 g to 25 kg.

Z		Low Pressure Chromato	graphy Reversed & Normal P	hases Characteristics
	Sorbent	Structure	Typical Characteristics	Typical Applications
	C18 PN: R33230B	S)~~~~~~~~	% C: ≥ 17.0 % Density: 0.639 g/mL	Purification of low to high polarity compounds     Reproducible purification without complexity and cost of preparative HPLC
	<b>C8</b> PN: R30830B	(S)~~~~	% C: ≥ 11.0 % Density: 0.586 g/mL	Less retention compared to C18     For highly hydrophobic pesticides, small peptides and large molecule drugs
9	Pentafluorophenyl (PFP) PN: R67530B	G	% C: ≥ 9.0 % Density: 0.761 g/mL	For an alternative selectivity, with aromatic ring interactions     For purification of conjugated compounds (isomers)
200	Phenyl (PHE) PN: R34030B	(s)—(¯)	% C: ≥ 8.0 % Density: 0.637 g/mL	Moderately non-polar sorbent     Alternative selectivity for aromatic compounds, compared to other reversed-phases
	<b>C4</b> PN: R32030B	3~~	% C: ≥ 8.0 % Density: 0.656 g/mL	Less retention compared to C18 and C8     For molecules with large hydrophobic regions
	<b>C1</b> PN: R33030B	Si-C <sub>1</sub>	% C: ≥ 5.0 % Density: 0.559 g/mL	Lower retention compared to other reversed-phases     For purification of polar and non-polar highly hydrophobic pharmaceutical products
	<b>Cyano</b> ( <i>CN</i> ) PN: R38030B	SI—≡N	% C: ≥ 7.0 % % N: ≥ 1.93 % Loading: ≥ 1.38 mmol/g Density: 0.703 g/mL	Versatile sorbent used either as normal or reversed-phas     Less polar than silica     For organic compounds with intermediate to extreme polarity
	Silica (Si) PN: R10030B	(si)—он	Density: 0.550 g/mL	Most popular sorbent for day-to-day use     For purification of non-ionic polar organic compounds
	Silica Premium PN: S10095D-A	Si — OH	Particle size: 25 μm Pore size: 90 Å Density: 0.450 g/mL	High performance sorbent for difficult separations (isomers     Higher loading capacity, faster flow rate, less solvent use
?	<b>Diol nec</b> PN: R35030B	(SI) OH OH	Loading: ≥ 0.97 mmol/g Density: 0.687 g/mL	For difficult separation of low to medium polarity samples     Can be used in HILIC mode     For mono and polysaccharides separation
200	<b>Amine</b> ( <i>NH</i> <sub>2</sub> , <i>WAX</i> ) PN: R52030B	Si NH <sub>2</sub>	Loading: ≥ 1.2 mmol/g Density: 0.700 g/mL	For purification of compounds with basic properties, or for monosaccharides separation
Normal	Acidic, Neutral & Basic Alumina PN: AUT-0053 AUT-0054 AUT-0055	$Al_2O_3$	Particle size: 75 - 150 μm Pore size: 70 Å	For aromatic compounds, aliphatic amines & compounds containing electronegative functions
	Florisil PN: AUT-0014	SiMgO <sub>3</sub>	Particle size: ≤ 75 µm Pore size: 80 Å	For separation of chlorinated pesticides, polychlorinated biphenyls ( <i>PCBs</i> ) & polysaccharides
	Silver Nitrate (AgNO <sub>3</sub> ) PN: R23530B	S + AgNO <sub>3</sub>	Loading: 10 % w/w Density: 0.604 g/mL	For separation of cis / trans isomers of unsaturated compounds (alkenes, lipids, steroids and terpenes)

If not otherwise stated, particle size is 40 - 63  $\mu m$  and pore diameter is 60 Å.

All phases are available endcapped and non-endcapped.

Other phases could be offered on a custom basis, contact us for more information: info@silicycle.com.






# **Typical Reversed and Normal Phases Applications**

The table below will help you select the right media to purify your compounds of interest. All phases are available either in bulk or pre-packed cartridges.

<u> </u>	Typical Applications Using Reversed and Normal Phases											
Analytes	Examples	C18	C8	C6	PFP	PHE	C4	C1	CN	NH <sub>2</sub>	Si	Diol
Biomolecules	Peptides, proteins	<b>√</b>	<b>√</b>	<b>✓</b>			<b>√</b>	<b>√</b>				<b>✓</b>
Nucleotides	Deoxyribonucleotides, ribonucleotides	<b>✓</b>								<b>✓</b>		
Lipids	Phospholipids		1	1			<b>1</b>	1		<b>✓</b>		
Carbohydrates	Sugars								1	<b>✓</b>		1
Glycosides	Glucosides, fructosides								1	1		1
Oligosaccharides	Malto-oligosaccharides									<b>√</b>		<b>√</b>
Pesticides	Organophosphates	1	1									
PCBs	Dichlorobiphenyl, trichlorobiphenyl	1			1	1						
PAHs	Anthracene, pyrene	1	1		1	1						
Drugs	Basic drugs, metabolites	1	1	1					1	1	1	
Alkaloids	Cocaine, morphine, nicotine, quinine	1	1						1		1	
Analgesics	Aspirin, acetaminophen, ibuprofen	1	1		1				1			
Cyclosporine	-	1							1			
Conjugated Compounds	Phenols, chloroanilines, steroids, caffeine	1	1	1	1	1	1	1				
Natural Compounds	Tannins, aflatoxins, flavonoids, carotenoids	1	1	1	1	1	1	1				
Fat-Soluble Vitamins	Vitamins A, D, E and K	1	1									
Water-Soluble Vitamins	Vitamins B and C									1	1	
Heterocyclic Compounds	Dioxins, furans	1										

AgNO<sub>3</sub> is particularly useful to separate isomers that present unsaturated groups. Neutral Alumina is used for the separation of aldehydes, ketones, quinines, esters, lactones and glucosides.

Florisil will help analyze pesticides, PCBs and PAHs.



# **Ion Exchange Phases Portfolio**

Y	<u> </u>	Low Pressi	ure Chromatography Ion I	Exchange Phases Characteristics
	Sorbent	Structure	Typical Characteristics	Typical Applications
	<b>Amine</b> ( <i>NH</i> <sub>2</sub> , <i>WAX</i> ) PN: R52030B	SI NH <sub>2</sub>	Loading: ≥ 1.2 mmol/g Density: 0.700 g/mL	<ul> <li>Weak anion exchanger (pK<sub>a</sub> of 9.8), positively charged at pH below 7.8</li> <li>For very strong anions (such as sulfonic acids), that may be too strongly retained on SAX phases</li> </ul>
	WAX-2 ( <i>Triethylamine</i> ) PN: R76530B	Si	Loading: ≥ 1.04 mmol/g Density: 0.761 g/mL	<ul> <li>Weak anion exchanger (pK<sub>a</sub> of 10.5), positively charged at pH below 8.5</li> <li>For catch &amp; release of compounds bearing a permanent negative charge (like salts of sulfonic acids)</li> </ul>
e Phases	SAX (TMA Chloride) PN: R66530B	SI N+ CI-	Loading: ≥ 0.90 meq/g Density: 0.700 g/mL	Strong anion exchanger, permanently positively charged (pH independant) For weak anions (such as carboxylic acids) that may not bind strongly enough on WAX phases For analysis of acidic drugs / analgesics, biomolecules & water-soluble vitamins
Exchange	SAX-2 nec (TMA Acetate) PN: R66430B	Si) N <sup>+</sup> CH <sub>3</sub> COO-	Loading: ≥ 0.71 mmol/g Density: 0.665 g/mL	Strong anion exchanger, with easily exchangeable acetate counter-ion (more than chloride ion)     For compounds with pK <sub>a</sub> < 5 (such as carboxylic acids)
lon	SCX (Tosic Acid) PN: R60530B	Si OH	Loading: ≥ 0.54 meq/g Density: 0.698 g/mL	• Strong cation exchangers (pK <sub>a</sub> < 1), permanently negatively charged (pH independant)
	SCX-2 ( <i>Propylsulfonic Acid</i> ) PN: R51230B	Si Si Si Si OH	Loading: ≥ 0.63 meq/g Density: 0.728 g/mL	For catch and release purification of weak cations (basic drugs / analgesics, biomolecules & water-soluble vitamins)
	WCX ( <i>Carboxylic Acid</i> ) PN: R70030B	Si	Loading: ≥ 0.92 mmol/g Density: 0.687 g/mL	<ul> <li>Weak cation exchanger (pK<sub>a</sub> of 4.8), neutralized at pH below 2.8</li> <li>For strong cationic species, that may bind too strongly on SCX phases</li> </ul>


For all sorbents, particle size is 40 - 63  $\mu m$  and pore diameter is 60 Å.

All bonded phases are available endcapped and non-endcapped.

Other phases could be offered on a custom basis, contact us for more information: info@silicycle.com.

### Counter-Ion Selectivity in Ion Exchange Mode

SAX phases are always paired with a counter-ion to neutralize the quaternary amine charge. But counter-ions have different selectivites and some are more easily removed from the silica gel by the analyte. You will find below the relative selectivity of standard counter-ions, compared to the hydroxyl ion OH<sup>-</sup> (lowest selectivity). Always choose a phase paired with a counter-ion less selective than the analyte.




17

Lower selectivity Easier to displace Harder to displace





# **Steps to Choose the Sorbent**



# SiliaBond Bulk Ordering Information

To build your own product number, just choose the right Codes for the desired Phase, Physical Characteristics, and Format: [Phase Code][Properties Code]-[Format Code]

Example: 100 g of C18 silica gel,  $\overline{40}$  - 63  $\mu$ m, 60  $\overline{\text{A}}$ : R33230B-100G.

SiliaBond phases are available on all irregular SiliaFlash silicas (R100-) and on all spherical SiliaSphere PC silicas (S100-). See page 11 for all available irregular & spherical silica types and corresponding codes.

You will find below the most common bare & bonded silica gels ordered in bulk.

Please note product numbers begin by R- for irregular silicas and by S- for spherical silicas.

Silia <i>Bond</i> P	hases
Phase	Code
Silia <i>Bond</i> Silica	R100
Silia <i>Bond</i> Amine	R520
Silia <i>Bond</i> Diol nec	R350
Silia <i>Bond</i> Cyano	R380
SiliaBond C18	R332
Silia <mark>Bond</mark> C8	R308
Silia <i>Bond</i> Phenyl	R340
Silia <i>Bond</i> PFP	R675
Silia <i>Bond</i> SCX	R605
SiliaBond SCX-2	R512
SiliaBond SAX nec	R665
SiliaBond SAX-2 nec	R664

SiliaBond Char	acteristics
Properties	Code
60 Å, 40 - 63 μm	30B
60 Å, 40 - 200 μm	40B
60 Å, 200 - 500 μm	70B
300 Å, 40 - 63 μm	30M
100 Å, 40 - 75 μm	30E
90 Å, 25 μm	95D-A

SiliaBond	Bulk Formats
Quantity	Code
5 g	5G
10 g	10G
25 g	25G
50 g	50G
100 g	100G
250 g	250G
500 g	500G
1 kg	1KG
5 kg	5KG
10 kg	10KG
25 kg	25KG

# **Flash Cartridges**

**Cartridge Design** 

Flash chromatography is one of the most used methods for compounds purification. Several types of flash chromatography are available, and it has been showed that the use of pre-packed flash cartridges improves purification efficiency compared to conventional flash, by offering superior reproducibility and productivity due to its tightly and homogeneously packed silica bed.

With SiliaSep, benefit from SiliCycle's renown quality: selectivity, speed and reliability.



#### Video: How flash chromatography works

Available in two silica gel grades, with various chemistries

- Spherical 25 µm (PREMIUM)
- Irregular 40 63 μm



19

### **Features and Benefits**

#### High silica gel quality, with low level of fines

- · No product contamination
- Homogeneous packing, no channelling (no peak tailing)
- High loading capacity (high surface area)
- · Direct transfer from TLC to flash chromatography

#### Reproducibility, reliability and safety

- Leak-free guaranteed by unique one-piece cartridge design
- Batch-to-batch reproducibility (stringent quality control)
- Excellent durability to withstand high pressures
- · Universal luer fittings for compatibility with any flash system

- Wide choice of cartridge sizes from 4 g to 3 Kg
- Purification scale-up from milligrams to kilograms
- Variety of sorbents to meet any separation needs

### Effective packing technology

- Consistent packing for reproducible high plate count (N)
- Excellent performance and separation
- High resolution with tight band definition (no tailing)
- · Great compound purity and recovery

#### Cost effectiveness

- Excellent performance / price ratio
- · Readily available, even for large volumes





# **Portfolio**

All our bare and bonded silica gels are available to be packed in SiliaSep flash cartridges to accomodate your chemistry.


	Silia <i>Sep</i> Flash Cartridges Adsorbents
Adsorbent Type	Adsorbent
Backbone	<ul> <li>Standard Silia<i>Flash</i> Irregular Silica, 40 - 63 μm, 60 Å</li> <li>PREMIUM Spherical Silica, 25 μm, 90 Å</li> <li>Acidic, Neutral and Basic Alumina, 50 - 75 μm, 55 Å</li> </ul>
Bonded phases	<ul> <li>SiliaBond Chromatographic Phases (reversed, normal and ion exchange phases)</li> <li>SiliaMetS Metal Scavengers (Thiol, DMT, etc.)</li> <li>SiliaBond Organic Scavengers (Amine, Tosic Acid, etc.)</li> </ul>



### **Formats**

\			Silias	ep Flash Cartridg	es Portfoli	0		
	Cartridge Format [Code]	Silica Weight (g)	Qty / Box	Dimensions (Diam. x Length*) (mm)	Column Volume ( <i>mL</i> )	Recommended Flow Rate (mL/min)	Loading Capacity (g)	Max Operating Pressure (psi / bar)
	SiliaSep 4 g [ISO04]	Bare: 4	20	12 x 98	4.9	15 - 25	Bare: 0.04 - 0.4	
R&D	3111a3cp 4 g [13004]	Bonded: ≥ 5	2	12 X 96	4.9	15 - 25	Bonded: 0.02 - 0.2	
and	SiliaSep 12 q [ISO12]	Bare: 12	20	21 x 117	17	20 - 40	Bare: 0.12 - 1.2	
Б >	3111436p 12 y [13012]	Bonded: ≥ 15	1	21 X 117	17	20 - 40	Bonded: 0.06 - 0.6	
Ver	SiliaSep 25 g [ISO25]	Bare: 25	15	21 x 165	31	20 - 45	Bare: 0.25 - 2.5	225 / 16
Discovery	3111 <b>436/ 23 y</b> [13023]	Bonded: ≥ 30	1	21 X 105	31	20 - 45	Bonded: 0.125 - 1.25	223 / 10
Dis	SiliaSep 40 g [ISO40]	Bare: 40	15	27 x 169	47	25 - 50	Bare: 0.4 - 4	
	3111836h 40 g [13040]	Bonded: ≥ 45	1	27 × 109	41	25 - 50	Bonded: 0.2 - 2	
	SiliaSep 80 g [ISO80]	Bare: 80	12	31 x 237	123	40 - 80	Bare: 0.8 - 8	
	Cinacop oo g [iccoo]	Bonded: ≥ 90	1	01 X 201	120	10 00	Bonded: 0.4 - 4	
	SiliaSep 120 g [IS120]	Bare: 120	10	36 x 256	190	60 - 120	Bare: 1.2 - 12	205 / 13
	omasop === g [.e==e]	Bonded: ≥ 130	1		200	00 120	Bonded: 0.6 - 6	200 / 20
SS	SiliaSep 220 q [IS220]	Bare: 220	4	60 x 195	306	60 - 190	Bare: 2.2 - 22	
əce		Bonded: ≥ 230	1				Bonded: 1.1 - 11	160 / 11
P.	SiliaSep 330 g [IS330]	Bare: 330	4	60 x 268	441	80 - 190	Bare: 3.3 - 33	100711
pu	Cinasch and a [10000]	Bonded: ≥ 360	1	00 X 200	771	00 100	Bonded: 1.65 - 16.5	
ıt a	SiliaSep XL 800 g** [IS750]	Bare: 800	2	78 x 382	1,500	200 - 300	Bare: 8 - 80	125 / 8
neı	31114367 AL 600 9 [13730]	Bonded: ≥ 870	1	76 X 362	1,500	200 - 300	Bonded: 4 - 40	12576
obi	SiliaSep XL 1,600 g** [11500]	Bare: 1,600	2	104 x 429	2,900	300 - 450	Bare: 16 - 160	100 / 7
Development and Process	31114364 VF 1,000 A [11200]	Bonded: ≥ 1,700	1	104 X 429	2,900	300 - 450	Bonded: 8 - 80	10077
De	SiliaSep XL 3 Kg** [ISO3KG]	Bare: 3,000	1	128 x 510		200 - 500	Bare: 30 - 300	100 / 7

<sup>\*</sup> Cartridge length includes luer-lock and connection tip.



# SiliaSep Flash Cartridges System Compatibility

Silia Sep cartridges are compatible with a variety of flash systems. They are a slip tip connection type, but they can be used with many leur lock systems. With some systems, and adapter kit might be required. The Table below is a guide to determine if an adapter is required to use SiliaSep cartridges with your flash system.

	SiliaSep System Compatibility	
System	SiliaSep Cartridges	Comments
Teledyne Isco™ CombiFlash®	COMPATIBLE	100 % compatible
Biotage Isolera™ & Selekt®	COMPATIBLE	100 % compatible
Büchi Pure™ & Sepacore™	COMPATIBLE	100 % compatible
Gilson PLC	COMPATIBLE	100 % compatible
Grace Reveleris™	COMPATIBLE	100 % compatible
Interchim PuriFlash™ & Spot II (Armen®)	COMPATIBLE	100 % compatible
Varian® (Analogix®) IntelliFlash® & SimpliFlash®	COMPATIBLE	100 % compatible
Biotage Horizon™	WITH ADAPTER	Use the Biotage Adapter Kit (PN: KAD-1006) of the Solvent Line Replacement (PN: KAD-1014
Biotage SP1 & SP4	WITH ADAPTER	Use Support Rings to allow the SiliaSep cartridge to sit on the instrument (Support Ring Kit PN: KAD-1008)
Biotage FlashMaster™	WITH ADAPTER	Use the FlashMaster Adapter Kit (PN: KAD- 1016) or connect a SiliaSep OT cartridge





20

Solutions for Chromatography and Purification

<sup>\*\*</sup> For SiliaSep XL formats, you may need to use an XL Adapter to connect the cartridge onto your system. Part number AUT-0127-2.

# **Method Development**


### Prediction of Column Volumes (CV)

TLC data can be used to predict flash purification, based on the relationship between TLC retention factor (Rf) and flash retention time (measured in column volume, CV). CV is the number of column volumes required to elute the component from the column, regardless of column dimensions.

So the first step to convert a TLC method in flash chromatography is to convert Rf into CV.

Rf and CV are inversely proportional: CV = 1 / Rf

The graph below shows that lower Rfs in TLC means greater CVs in flash (so better analyte retention). On the right is a chart giving CV values according to typical Rf values.





Video: Relationship between retention factor and column volume

As CV is a measure of analyte retention, then  $\Delta$ CV is a measure of two analytes separation and resolution:  $\Delta CV = CV_2 - CV_1 = (1 / Rf_2) - (1 / Rf_1)$ 

A	/						CV Val	ues Ac	cording	To Rf <sub>1</sub>	And Rf	<sub>2</sub> Value	s						
0.05	0.00																		
0.10	10.00	0.00																	
0.15	13.33	3.33	0.00																
0.20	15.00	5.00	1.67	0.00		,													
0.25	16.00	6.00	2.67	1.00	0.00														
0.30	16.67	6.67	3.34	1.67	0.67	0.00													
0.35	17.14	7.14	3.81	2.14	1.14	0.47	0.00												
0.40	17.50	7.50	4.17	2.50	1.50	0.83	0.36	0.00											
0.45	17.78	7.78	4.45	2.78	1.78	1.11	0.64	0.28	0.00										
0.50	18.00	8.00	4.67	3.00	2.00	1.33	0.86	0.50	0.22	0.00									
0.55	18.19	8.19	4.86	3.16	2.16	1.52	1.05	0.69	0.41	0.19	0.00								
0.60	18.35	8.35	5.02	3.35	2.35	1.68	1.21	0.85	0.57	0.35	0.16	0.00							
0.65	18.46	8.46	5.13	3.46	2.46	1.79	1.32	0.98	0.68	0.46	0.27	0.11	0.00						
0.70	18.60	8.60	5.27	3.60	2.60	1.93	1.46	1.10	0.82	0.60	0.41	0.25	0.14	0.00					
0.75	18.67	8.67	5.34	3.67	2.67	2.00	1.53	1.17	0.89	0.67	0.48	0.32	0.21	0.07	0.00				
0.80	18.75	8.75	5.42	3.75	2.75	2.08	1.61	1.25	0.97	0.75	0.56	0.40	0.29	0.15	0.08	0.00			
0.85	18.83	8.83	5.50	3.83	2.83	2.16	1.69	1.33	1.05	0.83	0.64	0.48	0.37	0.23	0.16	0.08	0.00		
0.90	18.90	8.90	5.57	3.90	2.90	2.23	1.76	1.40	1.12	0.90	0.71	0.55	0.44	0.30	0.23	0.15	0.07	0.00	
0.95	18.95	8.95	5.62	3.95	2.95	2.28	1.81	1.45	1.17	0.95	0.76	0.60	0.49	0.35	0.28	0.20	0.12	0.05	0.00
Rf <sub>1</sub> Rf <sub>2</sub>	0.05	0.10	0.15	0.20	0.25	0.30	0.35	0.40	0.45	0.50	0.55	0.60	0.65	0.70	0.75	0.80	0.85	0.90	0.95

# From TLC to Low Pressure Chromatography


It is now understood that TLC methods should be optimized so that compounds of interest elute with lower Rfs, ideally between 0.1 and 0.4, to maximize retention and analytes separation. To obtain these Rfs values, you can adjust the TLC solvent mixture by using different solvents with different polarities, and change the composition of the final TLC solvent mixture.

An optimized TLC method will assure you a better separation and purification of your compounds in low pressure chromatography, with optimal loading capacity (you will be able to load more on the cartridge if your compounds are well separated).

We recommend using a flash cartridge phase matching the TLC plate, for a more linear and easy method conversion. You should also run your flash chromatography with the same solvent conditions as your TLC method (in isocratic mode).

#### Case Study

We need to separate two analytes, 1 and 2. We will study two different TLC configurations.



#### **In summary**:

- The lower the Rfs, the greater ΔCV
- The greater the ΔCVs, the greater the separation and resolution between the spots (easier separation)

23


• The greater the  $\Delta$ CVs, the more sample can be loaded onto the column






# **Low Pressure Chromatography Loading Chart**

The chart below will help you choose the right cartridge size according to your sample size and your TLC results.



#### **Cartridge Size**



SILICYCLE (3)

The loading capacity depends on the sample itself, the column dimensions and the silica type. You will find below the sample loading we recommand with our SiliaSep flash cartridges. For easily separated compounds ( $\Delta CV > 6$ ) we suggest to load up to 5 % on irregular bonded phases, up to 10 % on bare irregular silica and up to 15 % on bare spherical silica.

$\Delta$			Low	Pressure	Chromato	graphy Lo	ading Cha	ırt				
							Loa	d (g)				
Dimensions	SiliaSep	SiliaSep	Diffic	cult Separa	ation	Medi	um Separ	ation '		Easy Se	paration	
ID x Length	Format	Phase	ΔCV = 0.1 - 0.6	ΔCV = 0.7 - 1.2	ΔCV = 1.3 - 1.8	ΔCV = 1.9 - 2.4	ΔCV = 2.5 - 3.1	ΔCV = 3.2 - 3.8	ΔCV = 3.9 - 4.5	ΔCV = 4.6 - 5.2	ΔCV = 5.3 - 6.0	ΔCV > 6
		Irregular Silica	0.040	0.080	0.120	0.160	0.200	0.240	0.280	0.320	0.360	0.400
12 x 98 mm	4 g	PREMIUM	0.052	0.104	0.156	0.208	0.260	0.312	0.364	0.416	0.468	0.520
		Bonded	0.020	0.040	0.060	0.080	0.100	0.120	0.140	0.160	0.180	0.200
		Irregular Silica	0.120	0.240	0.360	0.480	0.600	0.720	0.840	0.960	1.080	1.200
21 x 117 mm	12 g	PREMIUM	0.156	0.312	0.468	0.624	0.780	0.936	1.092	1.248	1.404	1.560
		Bonded	0.060	0.120	0.180	0.240	0.300	0.360	0.420	0.480	0.540	0.600
		Irregular Silica	0.250	0.500	0.750	1.000	1.250	1.500	1.750	2.000	2.250	2.500
21 x 165 mm	25 g	PREMIUM	0.325	0.650	0.975	1.300	1.625	1.950	2.275	2.600	2.925	3.250
		Bonded	0.125	0.250	0.375	0.500	0.625	0.750	0.875	1.000	1.125	1.250
		Irregular Silica	0.400	0.800	1.200	1.600	2.000	2.400	2.800	3.200	3.600	4.000
27 x 169 mm	40 g	PREMIUM	0.520	1.040	1.560	2.080	2.600	3.120	3.640	4.160	4.680	5.200
		Bonded	0.200	0.400	0.600	0.800	1.000	1.200	1.400	1.600	1.800	2.000
		Irregular Silica	0.800	1.600	2.400	3.200	4.000	4.800	5.600	6.400	7.200	8.000
31 x 237 mm	80 g	PREMIUM	1.040	2.080	3.120	4.160	5.200	6.240	7.280	8.320	9.360	10.400
		Bonded	0.400	0.800	1.200	1.600	2.000	2.400	2.800	3.200	3.600	4.000
		Irregular Silica	1.200	2.400	3.600	4.800	6.000	7.200	8.400	9.600	10.800	12.000
36 x 256 mm	120 g	PREMIUM	1.560	3.120	4.680	6.240	7.800	9.360	10.920	12.480	14.040	15.600
		Bonded	0.600	1.200	1.800	2.400	3.000	3.600	4.200	4.800	5.400	6.000
		Irregular Silica	2.200	4.400	6.600	8.800	11.000	13.200	15.400	17.600	19.800	22.000
60 x 195 mm	220 g	PREMIUM	2.860	5.720	8.580	11.440	14.300	17.160	20.020	22.880	25.740	28.600
		Bonded	1.100	2.200	3.300	4.400	5.500	6.600	7.700	8.800	9.900	11.000
		Irregular Silica	3.300	6.600	9.900	13.200	16.500	19.800	23.100	26.400	29.700	33.000
60 x 268 mm	330 g	PREMIUM	4.290	8.580	12.870	17.160	21.450	25.740	30.030	34.320	38.610	42.900
		Bonded	1.650	3.300	4.950	6.600	8.250	9.900	11.550	13.200	14.850	16.500
		Irregular Silica	8.000	16.000	24.000	32.000	40.000	48.000	56.000	64.000	72.000	80.000
78 x 382 mm	800 g	PREMIUM	10.400	20.800	31.200	41.600	52.000	62.400	72.800	83.200	93.600	104.000
		Bonded	4.000	8.000	12.000	16.000	20.000	24.000	28.000	32.000	36.000	40.000
		Irregular Silica	16.000	32.000	48.000	64.000	80.000	96.000	112.000	128.000	144.000	160.000
104 x 429 mm	1,600 g	PREMIUM	20.800	41.600	62.400	83.200	104.000	124.800	145.600	166.400	187.200	208.000
		Bonded	8.000	16.000	24.000	32.000	40.000	48.000	56.000	64.000	72.000	80.000
		1	Diffic	cult Separa	ation	Medi	um Separ	ation		Easy Se	paration	

For alumina sorbent, refer to the bare silica loading capacity.

It is worth noting that unlike with normal phases where it is easy to determine and optimize separation conditions using TLC plates, in reversed-phase chromatography the same cannot be said. First, the functions are impregnated on a TLC plate while they are bounded on SiliaBond products, thus the compounds do not have the exact same behaviors interacting with these products. Secondly, the plates are dry when the product is applied. However, the silica is already wetted in a chromatographic column. Since a reversed-phase is harder to wet than a normal phase, less contact area is available which tends to cause the products to migrate with the solvent front. Finally, the sample application is harder on a hydrophobic surface.

With this in mind, a general idea of the elution conditions can be found with some adjustments needed when transposing these onto a flash cartridge.



# SiliaSep Cartridges Cleaning and Re-Use

Pre-packed flash cartridges are designed and typically used for a single purification run (1-injection). Single-use gives the highest purification performance and the lowest solvent consumption. It is typically the easiest process to validate and it may give the lowest purification process cost.

It is possible to develop and validate a cleaning process that meets FDA requirements, so the flash cartridge can be used for multiple runs. This cleaning process is the client's responsibility. SiliCycle does not warranty any flash cartridge for multiple injections and all process validation is under the client's (sole) responsibility.

<b>\( \)</b>	Guidelines for Flash Cartridge Use in cGMP Environments
SiliaSep Phase	Recommended Use and Cleaning Procedure
Bare Silica (normal phase separations)	Porous silica is used in adsorption chromatography processes, where the product and its impurities "bind" to the surface with varying degrees of affinity. The solvent polarity is increased to desorb the product and its impurities at different elution volumes. While it is possible to elute nearly all the product from silica, some impurities typically remain at the end of each separation. If the cartridge is not fully cleaned, this remaining material may reduce the purification effectiveness and these impurities may elute in a subsequent run. Clearly, if the cartridge is planned to be used for a second or subsequent run, the process will require a validated cleaning protocol.  Some guidelines are given below:  • Single injection of a single batch of one API in this case, the cartridge is eluted and the purified product is collected. The cartridge is flushed and then discarded. This single-use process has the minimum solvent consumption and no-risk of cross-contamination.  • Multiple injections of a single batch of one API In this process, the full batch is too large to purify in a single run, and therefore multiple runs are required. Each injection is from a single batch or lot, and therefore the product and its impurities are identical in each injection or sample load. The cartridge must be cleaned between runs, but no cross-contamination is possible between batches.  Re-using silica cartridges for multiple injections within a single batch is a well accepted process decision. It must be demonstrated that each of the multiple injections gives the same elution profile and that the product purity is consistent in each of the sequential runs. Typically, the process control points are set to ensure that the impurity profile does not change more than 0.1 %.  This process can be modeled at the lab or pilot scale and then demonstrated at full production volume. In this process, the cleaning solvent is often 100 % of the most polar solvent in the elution mixture and is often carried out in reverse flow mode. The cartridg
C18 (reversed-phase separations)	Reversed-phase media is often used for multiple batches of a single API. However, due to the high cost and technical risk of batch-to-batch contamination, a fully validated cleaning procedure must be implemented. If a multiple lot strategy is nonetheless considered, the cleaning process will require a high level of data to support the decision.*  The cleaning protocol can be modeled at the lab or pilot scale and then demonstrated at full production volume. In this process, the cleaning solvent is often 100 % of the most polar solvent in the elution mixture ( <i>typically methanol</i> , <i>ethanol or acetonitrile</i> ), and is often carried out in reverse flow mode. The cartridge must be re-equilibrated, in normal flow mode, with the initial solvent conditions prior to the next injection. The cleaning step and re-equilibration step will each use a minimum of 3 column volumes of each solvent.
	It is recommended that C18 flash cartridges be dedicated to an individual API and never be used for multiple API compounds.

<sup>\*</sup> The data set must include analytical methods, such as HPLC and/or GC, and data to determine residue analysis. The standard assay is Total Organic Carbon (TOC) analysis. The upper and lower control limits must be set and defined for this process. The FDA does not set a number, but many organizations have used 1/1000 of a therapeutic dose of Product A in Product B as a guideline. This is a very challenging requirement, and the cost of cleaning solvents and time invested may be prohibitive.

26

# SILICYCLE (3)

# SiliaSep Cartridges Ordering Information

To build your own product number, just add the Format Code to the Phase PN: FLH-[Phase PN]-[Format Code] Example: SiliaSep C18 40 - 63 μm irregular grade, 4 g cartridge: FLH-R33230B-ISO04.

### SiliaSep Phases

### SiliaSep 40 - 63 µm Irregular Silica

Υ	SiliaSep P	hases	
Ph	ases	Phase PN	
	Bare Irregular Silica	R10030B	
	Acidic Alumina	AUT-0053	
	Neutral Alumina	AUT-0054	
	Basic Alumina	AUT-0055	-16
	C18	R33230B	
<u>S</u>	C8	R30830B	
S SIL	Phenyl	R34030B	
JLA	PFP	R67530B	
ZEGI	Amine	R52030B	
본	Diol nec	R35030B	VY
IZEI	Cyano	R38030B	
NAL	SCX	R60530B	
FUNCTIONALIZED IRREGULAR SILICA	SCX-2	R51230B	
Ž	SAX nec	R66530B	
_	SAX-2 nec	R66430B	•


### SiliaSep PREMIUM 25 µm Spherical Silica

Υ	SiliaSep PREI	MIUM Phases	
Pŀ	nases	Phase PN	
Т	Bare Spherical Silica	10095D-A	
	C18	03295D-A	
SILICA	C8	30895D-A	
	Phenyl	34095D-A	
ICAI	PFP	67595D-A	
SPHERICAL	Amine	52095D-A	
	Diol nec	35095D-A	
.IZE	Cyano	38095D-A	
NAL	SCX	60595D-A	
CTIC	SCX-2	51295D-A	
FUNCTIONALIZED	SAX nec	66595D-A	
_	SAX-2 nec	66495D-A	

Note: Other phases could be offered, like metal scavengers. Contact us for more information: info@silicycle.com.

### SiliaSep Formats

Y.		Silia <i>Sep</i> Forma	ts	
Formats	Qty / Box Bare Silica	Qty / Box Functi Standard Box	onalized Silica* Big Box	Format Code
4 g	20	2	10	ISO04
12 g	20	1	10	ISO12
25 g	15	1	10	ISO25
40 g	15	1	5	ISO40
80 g	12	1	5	ISO80
120 g	10	1	5	IS120
220 g	4	1	4	IS220
330 g	4	1	4	IS330
XL 800 g	2	1	-	IS750
XL 1,600 g	2	1	-	11500
XL 3 Kg	1	-	-	ISO3KG



· For bigger columns, please contact us.

Solutions for Chromatography and Purification

• For SiliaSep XL formats, you may need to use an XL Adapter to connect the cartridge onto your system. Part number: AUT-0127-2.







# SiliaSep Solid-Load Cartridges and Plungers

The use of solid-load technique (also called dry-load) is known to improve chromatography resolution, especially for compounds soluble only in strong solvents or in large volumes. Our plungers and solid-load cartridges are designed to allow injections in a dry form on SiliaSep flash cartridges.

### SiliaSep Solid-Load Cartridges

Four formats of solid-load cartridges are available: 10, 55, 150 and 700 mL (SiliCycle offers adapted plungers for the smaller ones - 10 and 55 mL- but not for the larger ones - 150 and 700 mL\*). The open top of the cartridge receives the bottom end of the plunger, and the luer-lock bottom connects to the flash cartridge.

Our solid-load cartridges are available either empty (to be filled with your own media) or pre-packed (with 2, 5, 10, 20, 65\* or 270\* g of *media*), whichever best suits your needs.

28

\* These two solid-load cartridges formats (150 and 700 mL) are compatible with the large-scale plunger of the CombiFlash Torrent® system.

### SiliaSep Solid-Load Cartridge Installation

#### A. EMPTY solid-load cartridges (to concentrate the sample and eliminate any solvent effect on the purification)

- 1. Mix your sample with bulk silica to make a silica-sample slurry. For a dry sample use a 1:1 ratio (1 g of silica for 1 g of dry sample) and for an oily sample use a 3:1 ratio (3 g of silica for 1 g of oily sample).
- 2. Dry the slurry by evaporating the solvent.
- Transfer the dried silica-sample powder in the empty solid-load cartridge and insert the lose frit.
- Connect the solid-load cartridge to the plunger (top) and the flash cartridge (bottom).



10 mL empty solid-load cartridge

#### B. PRE-PACKED solid-load cartridges (to allow direct injection of liquid samples)

- 1. Inject your liquid sample into the solid-load cartridge. You should be able to dilute your sample in 1 column volume at the most. If not, choose a bigger pre-packed solid-load cartridge.
- 2. Remove the dissolution solvent by gravity and / or heating.
- 3. Connect the solid-load cartridge to the plunger (top) and the flash cartridge (bottom).

These pre-packed solid-load cartridges are available with various SiliCycle's standard stationary phases, mainly: Silica, C18, Amine, Cyano and Diol (see Table on the following page for more information).



55 mL pre-packed solid-load cartridge




10 mL pre-packed solid-load cartridge

# SiliaSep Plungers

SiliCycle offers two plunger sizes: 10 and 55 mL. Connect the female luer-lock connection fitting to the flash system (solvent inlet) and insert the bottom part of the plunger inside the solid-load cartridge.

Both the compression plate and the perfluoroelastomer O-ring are solvent resistant, which permits the use of any mobile phase for the purification.



# SiliaSep Solid-Load Cartridges and Plungers Ordering Information

<u> </u>	SiliaSep Solid-Load Cartridges			
Product Number	Sorbent	Weight / Volume	Description	Qty / Box
AUT-0060-010	-	-	Plunger for 10 mL Solid Load Cartridge (16 mm)	1
AUT-0060-055	-	-	Plunger for 55 mL Solid Load Cartridge (23 mm)	1
SPL-R10030B-10U	Silica (40 - 63 μm)	2 g / 10 mL	SiliaSep Silica Pre-packed Solid-Load Cartridge, 2 g, 10 mL	20
SPL-R10030B-10X	Silica (40 - 63 μm)	5 g / 10 mL	SiliaSep Silica Pre-packed Solid-Load Cartridge, 5 g, 10 mL	20
SPL-R10030B-55Y	Silica (40 - 63 μm)	10 g / 55 mL	SiliaSep Silica Pre-packed Solid-Load Cartridge, 10 g, 55 mL	16
SPL-R10030B-55Z	Silica (40 - 63 μm)	20 g / 55 mL	SiliaSep Silica Pre-packed Solid-Load Cartridge, 20 g, 55 mL	16
SPL-R10030B-065*	Silica (40 - 63 μm)	65 g / 150 mL	SiliaSep Silica Pre-packed XL Solid-Load Cartridge, 65 g, 150 mL*	12
SPL-R10030B-270*	Silica (40 - 63 μm)	270 g / 700 mL	SiliaSep Silica Pre-packed XL Solid-Load Cartridge, 270 g, 700 mL*	6
SPL-R52030B-10X	Amine	5 g / 10 mL	SiliaSep Amine Pre-packed Solid-Load Cartridge, 5 g, 10 mL	20
SPL-R52030B-55Z	Amine	20 g / 55 mL	SiliaSep Amine Pre-packed Solid-Load Cartridge, 20 g, 55 mL	16
SPL-R35030B-10X	Diol	5 g / 10 mL	SiliaSep Diol Pre-packed Solid-Load Cartridge, 5 g, 10 mL	20
SPL-R35030B-55Z	Diol	20 g / 55 mL	SiliaSep Diol Pre-packed Solid-Load Cartridge, 20 g, 55 mL	16
SPL-R38030B-10X	Cyano	5 g / 10 mL	SiliaSep Cyano Pre-packed Solid-Load Cartridge, 5 g, 10 mL	20
SPL-R38030B-55Z	Cyano	20 g / 55 mL	SiliaSep Cyano Pre-packed Solid-Load Cartridge, 20 g, 55 mL	16
SPL-R33230B-10X	C18	5 g / 10 mL	SiliaSep C18 Pre-packed Solid-Load Cartridge, 5 g, 10 mL	20
SPL-R33230B-55Z	C18	20 g / 55 mL	SiliaSep C18 Pre-packed Solid-Load Cartridge, 20 g, 55 mL	16
SPL-0009-010	Empty	- / 10 mL	SiliaSep Empty Solid-Load Cartridge, 10 mL (with 200 frits)	100
SPL-0012-055	Empty	- / 55 mL	SiliaSep Empty Solid-Load Cartridge, 55 mL (with 200 frits)	100
AUT-0090-150*	Empty	- / 150 mL	SiliaSep Empty Solid-Load Cartridge, 150 mL (with 24 frits)*	12
AUT-0090-700*	Empty	- / 700 mL	SiliaSep Empty Solid-Load Cartridge, 700 mL (with 12 frits)*	6

<sup>\*</sup> These two solid-load cartridges formats (150 and 700 mL) are compatible with the large-scale plunger of the CombiFlash Torrent® system.

Notes: For optimal purification performance, solvent removal under vacuum is highly recommended.

Solutions for Chromatography and Purification

Other phases can be offered pre-packed in our solid-load cartridges, contact us for more information: info@silicycle.com





# SiliaSep OT (Open Top) Flash Cartridges

SiliaSep OT cartridges are mainly used with vacuum manifolds and automated SPE equipments. They are also directly compatible with FlashMaster™ systems.

### **Ordering Information**

To build your own product number, just add the Format Code to the Phase PN: FLH-[Phase PN]-[Format Code] Example: SiliaSep OT C18, 15 g cartridge: FLH-R00230B-70i.

SiliaSe	OT Phases
Phase	Phase PN
Bare Silica	R10030B
C18	R00230B
C8	R31030B
Phenyl	R34030B
PFP	R67530B
Amine	R52030B
Diol nec	R35030B
Cyano	R38030B
SCX	R60530B
SCX-2	R51230B
SAX nec	R66530B
SAX-2 nec	R66430B

SiliaSep OT Formats				
Format	Qty / Box	Dimensions (ID x length) (mm)	Forma Bare Silica	at Code Functionalized Silica
2 g / 12 mL	20	15.8 x 90	<b>15</b> U	SPE-[Phase PN]-12U
5 g / 25 mL	20	20.5 x 100	25X	SPE-[Phase PN]-20X
10 g / 70 mL	16	26.8 x 154	70Y	
15 g / 70 mL	16	26.8 x 154	70i	
20 g / 70 mL	16	26.8 x 154	70Z	
25 g / 150 mL	10	38.2 x 170	95K	
50 g / 150 mL	10	38.2 x 170	95M	
70 g / 150 mL	10	38.2 x 170	95N	

#### Notes

- Other phases could be offered, contact us for more information: info@silicycle.com.
- SiliaSep OT cartridges are also available with bar code for automation purposes
- Maximum operating pressure: 60 psi.



# Thin Layer Chromatography (TLC)

- Wide choice of sizes, sorbents, and thicknesses available
- Excellent reproducibility between Silia *Plate* TLC plates and bulk silicas or flash cartridges

The hardness of our silica layer, combined to a homogeneous coating and layer thickness, allows excellent separations. Each TLC batch is chemically and physically controlled by our Quality Control department to ensure lot-to-lot and layer-to-layer reproducibility.

# SiliaPlate TLC Plates

Thin-layer chromatography (TLC) is a quick, simple and inexpensive analytical technique frequently used in various laboratories as it is one of the most verstatile. It is used for reaction monitoring, screening, and compound purity evaluation.

Rapid and cost-efficient selection and optimization of chromatographic conditions prior to flash chromatography purification or HPLC analysis.

Besides speed and low cost, TLC analysis presents other non-negligible advantages like the small quantity of compound required and high sample throughput capability (up to 20 samples simultaneously).

Like column chromatography, TLC is a solid-liquid partitioning technique, in which the sample is applied to the plate as a small spot near the base of the plate. The moving liquid phase is then allowed to ascend the plate, causing the sample to partition between moving and stationary phase.

### SiliaPlate Features and Benefits

For more than 20 years, SiliCycle has been offering a wide selection of TLC plates in various sizes (plate size, thickness, backing) and chemistries (10 % Silver Nitrate, CN, C18, NH<sub>2</sub>). SiliaPlate represents an efficient and economical alternative to other TLC plate manufacturers while demonstrating high separation power, which is due to our narrow particle size distribution silica gel.

The extraordinary silica layer hardness combined to a homogeneous coating and layer thickness allows excellent separation. Each TLC batch is chemically and physically controlled by our Quality Control department to ensure lot-to-lot and layer-to-layer reproducibility.

# **Selection Guide**

# **Plate Types**

SiliCycle offers different types of plates for thin-layer chromatography applications: classical TLC, high performance TLC (also called HPTLC) and preparative TLC (PLC). The plate types are selected based on the analysis required and the available budget.

Differences Between Classical TLC, HPTLC and PLC				
Properties	Classical TLC	HPTLC	Preparative (PLC)	
Applications	Quick, inexpensive, flexible and classical separations	Highly sophisticated separations, complex samples	Purification on a TLC plate	
Analysis	Qualitative	Qualitative & quantitative	Quantitative	
Detection	UV - Stains	Instrumented analysis (use of scanners for detection)	UV	
Distribution [Mean Particle Size]	5 - 20 μm [ <i>10 - 14 μm</i> ]	4 - 8 μm [5 - <i>6 μm</i> ]	5 - 40 μm [22 - 25 μm]	
Layer Thickness	200 - 250 μm	150 - 200 μm	500 - 2,000 μm	
Typical Sample Volume	1 - 5 μL	0.1 - 0.5 μL	5 - 20 μL	

31



Solutions for Chromatography and Purification





### **TLC Backings**

TLC plates are available with different backings: rigid (glass-backed) and flexible (aluminum & plastic-backed).

Y.	TLC Backings Comparison				
Properties	Glass	Aluminum	Plastic		
Advantages	Rigid     High chemical resistance     High heating stability and charring resistance     Transparent	Thin and low fragility Low weight & consequent shipping costs High heating stability Possible to cut with scissors Can be stored in notebook	Thin Low fragility Possible to cut with scissors High chemical resistance Can be stored in notebook		
Disadvantages	Thick and high fragility Impossible to cut with scissors Cannot be stored in lab notebook High weight & consequent shipping costs Large shelf space	Low chemical resistance     Opaque	Medium weight     Opaque     Heating stability up tp 175°C     Possible cracking of matrix due to high flexibility		
Thickness (approx.)	2.0 - 2.5 mm	1.5 - 2.0 mm	1.5 - 2.0 mm		
Heating Stability	High	High	Below 175°C		
Fragility	High	Low	Low		
Cutting with Scissors	Impossible	Easily	Possible		
Chemical Resistance	High	Low	High		

### **Layer Thicknesses**

The layer thickness is related to the nature of the analysis (analytical or preparative) as well as the performance of the plate (TLC or HPLTC). The most common layer thicknesses are 150 - 200 µm (HPTLC plates), 200 - 250 µm (analytical TLC plates), and 500 - 2,000 µm (preparative TLC plates)

### **UV** Indicator

Silia Plate TLC Plates could contain inorganic fluorescent indicator for UV detection of colorless substances: F254 (manganese activated zinc silicate with absorption maximum at 254 nm; green fluorescence). Therefore, samples which absorb short-wave UV at 254 nm (F254) can be viewed under UV light. As a matter of fact, substances absorbing in the respective region of wavelength cause a quenching of the fluorescence and appear as dark spots on the TLC.

#### Channeled

Those TLC plates have 'channels' (alternating zones with and without silica) which allows to separate several samples and prevent cross contamination from a sample to another.

### **Preadsorption Zone**

Silia Plate TLC Plates with preadsorption zone allow to load quickly and easily samples, even large volumes of diluted samples. The adsorbent in the preadsorption zone is a large pore concentrating adsorbent and the other one is a selective layer for separation.

With this kind of plates, the sample always concentrate in a narrow band (no matter the way sample is loaded; spots do not have to be the exact same size or on the exact same line).

#### Binder

SiliCycle offers two types of binder, with different sensitivities and areas of applications: 'B' and 'BK'

- B's layer is polymeric: a small percentage of inorganic, hardening agent has been added for a uniform and hard surface, smooth and dense, that will not crack, blister or swell up. They were designed for maximum robustness of the binder: they are very easy to handle and to write on, as well as completely wettable. They are compatible with all solvents, yet, they might oxidize a bit faster when dipped into KMnO, (fading in a few minutes from flashy purple to yellow ocher). Also, spots are a bit less definite when using CAM as a revelatory. Such binder also contains a higher percentage of fluorescent indicator for greater brilliance of spots and less background noise from the silica layer.
- <u>BK</u>'s layer is gypsum (calcium sulfate), and do not contain the polymeric additive that provides the former plates a harder surface and ruggedness. This means that the layer is softer, so spots can be easily scrapped off from the glass support, and are particularly recommended for aggressive visualization methods (strong charring, CAM staining solution) or, if dipped into KMnO<sub>4</sub>, ought to remain bright-purple a longer period of time.

Here is a chart which can hopefully help you quickly choose the right plate for your specific application.

SiliaPlate TLC Plates Binders				
TLC Plate Binder	В	вк		
Example	TLG-R10014B-323	TLG-R10014 <mark>BK-</mark> 323		
UV Fluorescence (F <sub>254</sub> )	Higher brightness Less background noise from layer	Yes		
	Stable in almost a	all organic solvents		
Binder Sensitivity	Increased separation efficiency	Resistant to aggressive visualization methods		
Surface Layer	Robust and rugged	Easily scratched off		
Water Tolerance	Up to 80 %	Up to 40 %		
Specific Surface (BET)	≈ 500	) m²/g		
Mean Pore Size	60 Å			
Mean Pore Volume	0.75	mL/g		
Distribution (Mean Particle Size)	5 - 20 μm [	10 - 14 μm]		
Layer Thickness	≈ 25	50 μm		
Stain Compatibility				
KMnO <sub>4</sub>	Compatible	Highly compatible		
CAM	Сотр	patible		
p-Anisaldehyde	Compatible	Highly compatible		
Ninhydrin	Highly co	ompatible		
Vanilin	Highly co	ompatible		

### **Sorbents**

Classical Silica Gel	Reversed-Phases	Special Phases
A universal matrix for daily, fast, reliable analysis of the largest spectra of molecules  The particle size distribution used for the silica is related to the nature of the plate.	'' '	s polar than the stationary phase. Inversely, in reversed
For standard TLC, silica gel with a mean particle size of 10 - 14 $\mu m$ is used compared to HPTLC, where a smaller particle size is required. In both cases, pore diameter is always 60 Å.	When satisfactory separations cannot be achieved b been designed for specific applications:	y unmodified silica, other functionalized matrices have
	C2, C8 and C18 phases are functionalizations of silica performed using organosilanes of various chain lengths. Retention of molecules & ability to tolerate water in the moving phase are directly dependent on the chains length.	Diol and Cyano ( <i>CN</i> ) are moderately polar. They can thus be suitable for both normal and reversed-phase chromatography, depending on your application.  Amino phases ( <i>NH</i> <sub>2</sub> ) show weak anion exchange characteristics, great for charged compounds.

**Available Sorbents** 



### **Matrices** (or Adsorbents)

Various adsorbents can be used for TLC coating; silica, aluminum oxide, florisil, etc. For more than 20 years, SiliCycle has been offering a wide selection of TLC plates in various sizes (plate size, thickness, backing) and chemistries (Silver Nitrate, CN, C18, NH). More than 80 % of all purifications are performed using silica gel as the adsorbent.

Available Matrices			
Silica Gel	Aluminum Oxide		
Can be unmodified or functionalized, and suitable for a myriad of molecules of functionalities & polarities, such as aflatoxins, alkaloids, barbiturates, fatty acids, flavonoids, glycosides, lipids, nucleosides, proteins, pesticides, sweeteners, vitamins and so on.	Aluminum oxide (commonly called Alumina) is the second most commonly used matrix, and shows similar selectivity to that of silica, with 3 different pH ranges (basic, neutral, acidic).  Popular applications include the separation for alkaloids, aliphatic compounds, aromatics, steroids, etc.		

# **Plate Sizes**

4	Available Sizes				
Standard TLC Plates	Micro TLC Plates	Scored TLC Plates			
Silia Plate TLC plates are available in the following standard sizes depending on the coating used:  20 x 20 cm  10 x 20 cm  5 x 20 cm  5 x 10 cm  Example:	Also for your convenience, SiliCycle provides ready-to-use micro TLC plates in the following formats:  2.5 x 10 cm  2.5 x 7.5 cm  Example:	An interesting compromise between standard and micro plate sizes is our Scored Silia Plate (glass backing).  Three different formats are available and possible cut combinations are shown in the image below.  20 x 20 cm plates scored to four 5 x 20 cm plates (or multiple of 5 cm width)  10 x 20 cm plates scored to eight 2.5 x 10 cm plates (or multiple of 2.5 cm width)  5 x 20 cm plates scored to eight 2.5 x 5 cm plates (or multiple of 2.5 cm width)			
		25cm 25cm 25cm 25cm 25cm 25cm 25cm 25cm			

# Silia Plate Ordering Information

Please note that this is an overview of plates that SiliCycle offers.

Different sizes are available, as well as more exotic layers for special separations (chiral layers, layers for surfactant separations, for PAH analysis, layers for basic or acidic ion exchange, cellulose layers, etc.). Contact us: info@silicycle.com.

Various combinations are possible with Silia Plate TLC plates and are summarized in the table below.

SiliaPlate TLC Plates Portfolio				
Properties	Analytical	HPTLC	Preparative	
Available Backings				
Glass	Yes	Yes	Yes	
Aluminum	Yes	Yes	No	
Plastic	Yes	No	No	
Available Adsorbents				
Bare Silica	Yes	Yes	Yes	
Functionalized Silica	No	Yes	Yes	
Silica Specifications		-		
Mean Particle Size	10 - 14 μm	5 - 6 μm	22 - 25 μm	
Mean Pore Diameter	60 Å	60 Å	60 Å	
Type of Plate Available				
Scored Plate	Yes	No	Yes	
Channeled Plate	Yes	No	No	
Layer Thickness	Glass: 250 μm Flexible: 200 μm	Glass: 150 - 200 μm	Glass: • 500 μm • 1,000 μm  Flexible: • 1,500 μm • 2,000 μm	
Plate Size	<ul> <li>2.5 x 5 cm</li> <li>2.5 x 7.5 cm</li> <li>2.5 x 10 cm</li> <li>5 x 10 cm</li> <li>5 x 20 cm</li> <li>10 x 20 cm</li> <li>20 x 20 cm</li> </ul>	<ul> <li>2.5 x 5 cm</li> <li>2.5 x 7.5 cm</li> <li>2.5 x 10 cm</li> <li>5 x 10 cm</li> <li>5 x 20 cm</li> <li>10 x 20</li> <li>20 x 20 cm</li> </ul>	• 20 x 20 cm	



Video: Sample mass loading rule of thumb for 20 x 20 cm plates





# **Glass-backed TLC Plates**

	Glass-bac	ked Analytical TLC Plates			
PN	Sorbent	Plate Size	Thickness	UV Indicator	Qty / Bo
Silica					
TLG-R10014B-124	Silica, Hard Layer	2.5 x 7.5 cm	250 μm	F254	100
TLG-R10014B-424	Silica, Hard Layer	5 x 20 cm	250 μm	F254	100
TLG-R10014B-323	Silica, Hard Layer	20 x 20 cm	250 μm	F254	25
TLG-R10014B-323N	Silica, Hard Layer	20 x 20 cm	250 μm	None	25
TLG-R10014BK-417	Silica, optimized for KMnO <sub>4</sub> revelation	2.5 x 5 cm	250 μm	F254	200
TLG-R10014BK-124	Silica, optimized for KMnO <sub>4</sub> revelation	2.5 x 7.5 cm	250 μm	F254	100
TLG-R10014BK-527	Silica, optimized for KMnO <sub>4</sub> revelation	5 x 10 cm	250 μm	F254	200
TLG-R10014BK-424	Silica, optimized for KMnO <sub>4</sub> revelation	5 x 20 cm	250 μm	F254	100
TLG-R10014BK-725	Silica, optimized for KMnO <sub>4</sub> revelation	10 x 20 cm	250 μm	F254	50
TLG-R10014BK-323	Silica, optimized for KMnO <sub>4</sub> revelation	20 x 20 cm	250 μm	F254	25
TLG-R10014BK-323N	Silica, optimized for KMnO <sub>4</sub> revelation	20 x 20 cm	250 μm	None	25
TLG-R10014BKB-323	Silica, optimized for KMnO <sub>4</sub> revelation	20 x 20 cm	250 μm	F254 , F366	25
Channeled with preadso			•	'	
TLGCZ-R10011B-423	Silica	5 x 20 cm	250 μm	F254	25
TLGCZ-R10011B-723	Silica	10 x 20 cm	250 μm	F254	25
TLGCZ-R10011B-323	Silica	20 x 20 cm	250 μm	F254	25
TLGCZ-R10011B-323N	Silica	20 x 20 cm	250 μm	None	25
Scored TLC plates		'			
TLGSR10011B-423	Silica	5 x 20 cm, scored to 2.5 x 5 cm	250 μm	F254	25
TLGSR10011B-424	Silica	5 x 20 cm, scored to 2.5 x 5 cm	250 μm	F254	100
TLGSR10011B-723	Silica	10 x 20 cm, scored to 2.5 x 10 cm	250 μm	F254	25
TLGSR10011B-323	Silica	20 x 20 cm, scored to 5 x 20 cm	250 μm	F254	25
Functionalized silica & c	other adsorbents	'	•	•	
TLG-AUT0014-423	Florisil	5 x 20 cm	250 μm	F254	25
TLG-AUT0014-723	Florisil	10 x 20 cm	250 μm	F254	25
TLG-AUT0014-323	Florisil	20 x 20 cm	250 μm	F254	25
TLG-AUT0337-323B	Basic Alumina	20 x 20 cm	250 μm	F254	25
TLG-AUT0337B-424N	Neutral Alumina	5 x 20 cm	250 μm	F254	100
TLG-AUT0337-323N	Neutral Alumina	20 x 20 cm	250 μm	F254	25
TLG-AUT0337-323NF	Neutral Alumina	20 x 20 cm	250 μm	None	25
TLG-AUT0337B-323N	Neutral Alumina	20 x 20 cm	250 μm	F254	25







# **Glass-backed TLC and HPTLC Plates**


	Glass-backed Anal	ytical TLC and HPTLC Plates			
PN	Sorbent	Plate Size	Thickness	UV Indicator	Qty / Box
Functionalized silica &	other adsorbents				
TLG-R30314BK-213	C18 (100 %), optimized for KMnO <sub>4</sub> revelation	10 x 10 cm	150 μm	F <sub>254</sub>	25
TLG-R30314BK-213N	C18 (100 %), optimized for KMnO <sub>4</sub> revelation	10 x 10 cm	150 μm	None	25
TLG-R30411B-213	C18 (13 %)	10 x 10 cm	150 μm	F <sub>254</sub>	25
TLG-R30411B-303	C18 (13 %)	20 x 20 cm	150 μm	F <sub>254</sub>	25
TLG-R30414B-313	C18 (13 %)	20 x 20 cm	200 μm	F <sub>254</sub>	25
TLG-R30411B-323	C18 (13 %)	20 x 20 cm	250 μm	F <sub>254</sub>	25
TLG-R31011B-203	C8	10 x 10 cm	150 μm	F <sub>254</sub>	25
TLG-R31011B-303	C8	20 x 20 cm	150 μm	F <sub>254</sub>	25
TLG-R32611B-203	C2	10 x 10 cm	150 μm	F <sub>254</sub>	25
TLG-R32611B-303	C2	20 x 20 cm	150 μm	F <sub>254</sub>	25
TLG-R32614BK-313	C2, optimized for KMnO <sub>4</sub> revelation	20 x 20 cm	200 μm	F <sub>254</sub>	25
TLG-R32614BK-713	C2, optimized for KMnO <sub>4</sub> revelation	10 x 20 cm	200 μm	F <sub>254</sub>	25
TLG-R35011B-713	Diol	10 x 20 cm	200 μm	F <sub>254</sub>	25
TLG-R35014BK-213	Diol, optimized for KMnO <sub>4</sub> revelation	10 x 10 cm	200 μm	F <sub>254</sub>	25
TLG-R35014BK-713	Diol, optimized for KMnO <sub>4</sub> revelation	10 x 20 cm	200 μm	F <sub>254</sub>	25
TLG-R35014BK-313	Diol, optimized for KMnO <sub>4</sub> revelation	20 x 20 cm	200 μm	F <sub>254</sub>	25
TLG-R38011B-203	Cyano (CN)	10 x 10 cm	150 μm	F <sub>254</sub>	25
TLG-R38011B-723	Cyano (CN)	10 x 20 cm	150 μm	F <sub>254</sub>	25
TLG-R38011B-303	Cyano (CN)	20 x 20 cm	150 μm	F <sub>254</sub>	25
TLG-R38014BK-213	Cyano (CN), optimized for KMnO <sub>4</sub> revelation	10 x 10 cm	200 μm	F <sub>254</sub>	25
TLG-R38014BK-713	Cyano (CN), optimized for KMnO <sub>4</sub> revelation	10 x 20 cm	200 μm	F <sub>254</sub>	25
TLG-R38014BK-313	Cyano (CN), optimized for KMnO <sub>4</sub> revelation	20 x 20 cm	200 μm	F <sub>254</sub>	25
TLG-R52011B-203	Amine (NH <sub>2</sub> )	10 x 10 cm	150 μm	F <sub>254</sub>	25
TLG-R52011B-723	Amine (NH <sub>2</sub> )	10 x 20 cm	150 μm	F <sub>254</sub>	25
TLG-R52011B-303	Amine (NH <sub>2</sub> )	20 x 20 cm	150 μm	F <sub>254</sub>	25
TLG-R52014BK-213	Amine (NH <sub>2</sub> ), optimized for KMnO <sub>4</sub> revelation	10 x 10 cm	200 μm	F <sub>254</sub>	25
TLG-R52014BK-713	Amine (NH <sub>2</sub> ), optimized for KMnO <sub>4</sub> revelation	10 x 20 cm	200 μm	F <sub>254</sub>	25
TLG-R52014BK-313	Amine (NH <sub>2</sub> ), optimized for KMnO <sub>4</sub> revelation	20 x 20 cm	200 μm	F <sub>254</sub>	25
TLG-R23511B-423	AgNO <sub>3</sub> (10 %)	5 x 20 cm	250 μm	F <sub>254</sub>	25
TLG-R23511B-303	AgNO <sub>3</sub> (10 %)	20 x 20 cm	250 μm	F <sub>254</sub>	25
TLG-R23611B-423	AgNO <sub>3</sub> (15 %)	5 x 20 cm	250 μm	F <sub>254</sub>	25
TLG-R23611B-323	AgNO <sub>3</sub> (15 %)	20 x 20 cm	250 μm	F <sub>254</sub>	25
TLG-R23711B-423	AgNO <sub>3</sub> (20 %)	5 x 20 cm	250 μm	F <sub>254</sub>	25
TLG-R23711B-323	AgNO <sub>3</sub> (20 %)	20 x 20 cm	250 μm	F <sub>254</sub>	25
TLG-R23M11B-323	AgNO <sub>3</sub> (5-10-15-20 %, 5 TLC each)	5 x 20 cm	250 μm	F <sub>254</sub>	5 x 4
TLGSR1234511B-723	Trial Packing of Functionalized Silica	10 x 20 cm, scored to 2.5 x 10 cm	150 μm	F <sub>254</sub>	25
TLG-AUT0308-203	RP Silanized	10 x 10 cm	150 μm	F <sub>254</sub>	25





# **Glass-backed Preparative TLC Plates**

	Glass-backed	Preparative TLC Plates			
PN	Sorbent	Plate Size	Thickness	UV Indicator	Qty / Box
Silica					1
TLG-R10011B-333	Silica	20 x 20 cm	500 μm	F <sub>254</sub>	25
TLG-R10011B-341	Silica	20 x 20 cm	1,000 μm	F <sub>254</sub>	25
TLG-R10011B-363	Silica	20 x 20 cm	1,500 μm	F <sub>254</sub>	25
TLG-R10011B-353	Silica	20 x 20 cm	2,000 μm	F <sub>254</sub>	25
Functionalized silica & o	ther adsorbents			-	
TLG-AUT0337-343N	Neutral Alumina	20 x 20 cm	1,000 μm	F <sub>254</sub>	25
TLG-AUT0337-343NF	Neutral Alumina	20 x 20 cm	1,000 μm	None	25
TLG-AUT0337-443	Neutral Alumina	5 x 20 cm	1,000 μm	F <sub>254</sub>	25
TLG-AUT0337-443F	Neutral Alumina	5 x 20 cm	1,000 μm	None	25
TLG-AUT0337B-341N	Neutral Alumina	20 x 20 cm	1,000 μm	None	15
TLG-R23511B-433	AgNO <sub>3</sub> (10 %)	5 x 20 cm	500 μm	F <sub>254</sub>	25
TLG-R23511B-333	AgNO <sub>3</sub> (10 %)	20 x 20 cm	500 μm	F <sub>254</sub>	25
TLG-R23611B-433	AgNO <sub>3</sub> (15 %)	5 x 20 cm	500 μm	F <sub>254</sub>	25
TLG-R23611B-333	AgNO <sub>3</sub> (15 %)	20 x 20 cm	500 μm	F <sub>254</sub>	25
TLG-R23711B-433	AgNO <sub>3</sub> (20 %)	5 x 20 cm	500 μm	F <sub>254</sub>	25
TLG-R23711B-333	AgNO <sub>3</sub> (20 %)	20 x 20 cm	500 μm	F <sub>254</sub>	25
TLG-R30411B-341	C18 (13 %)	20 x 20 cm	1,000 μm	F <sub>254</sub>	15
TLG-R30414BK-341	C18 (15 %), optimized for KMnO <sub>4</sub> revelation	20 x 20 cm	1,000 μm	F <sub>254</sub>	15
Scored preparative TLC	plates				
TLGSR10011B-333	Silica	20 x 20 cm, scored to 5 x 20 cm	500 μm	F <sub>254</sub>	25
TLGSR10011B-341	Silica	20 x 20 cm, scored to 5 x 20 cm	1,000 μm	F <sub>254</sub>	25
TLGSR10011B-353	Silica	20 x 20 cm, scored to 5 x 20 cm	2,000 μm	F <sub>254</sub>	25



# **Aluminum-backed TLC Plates**

	Aluminum-back	ed Analytical TLC Plates			
PN	Sorbent	Plate Size	Thickness	UV Indicator	Qty / Box
Silica					Î
TLA-R10011B-005	Silica	4 x 8 cm	150 μm	F <sub>254</sub>	50
TLA-R10011B-124	Silica	2.5 x 7.5 cm	200 μm	F <sub>254</sub>	200
TLA-R10011B-323	Silica	20 x 20 cm	200 μm	F <sub>254</sub>	25
TLA-R10011B-323N	Silica	20 x 20 cm	200 μm	None	25
TLA-R10011B-415	Silica	5 x 20 cm	200 μm	F <sub>254</sub>	50
TLA-R10011B-515	Silica	5 x 10 cm	200 μm	F <sub>254</sub>	50
TLA-R10011B-712	Silica	10 x 20 cm	200 μm	F <sub>254</sub>	20
TLA-R10014BK-1112	Silica, optimized for KMnO <sub>4</sub> revelation	5 x 7.5 cm	200 μm	F <sub>254</sub>	20
Functionalized silica & o	ther adsorbents				
TLA-AUT0337-323N	Neutral Alumina	20 x 20 cm	200 μm	F254	25
TLA-AUT0337-323NF	Neutral Alumina	20 x 20 cm	200 μm	None	25
TLA-R30411B-005	Silica C18 (13 %)	4 x 8 cm	150 μm	F254	50
TLA-R30411B-303	Silica C18 (13 %)	20 x 20 cm	150 μm	F254	25
TLA-R30411B-405	Silica C18 (13 %)	5 x 20 cm	150 μm	F254	50
TLA-R30414BK-303	Silica C18 (13 %), opt. for KMnO <sub>4</sub> revelation	20 x 20 cm	150 μm	F254	25
TLA-R52014BK-005	Amine $(NH_2)$ , optimized for $KMnO_4$ revelation	4 x 8 cm	150 μm	F254	50







### **Plastic-backed TLC Plates**

	Plastic-backed Analytical TLC Plates					
PN	Sorbent	Plate Size	Thickness	UV Indicator	Qty / Box	
Silica						
TLP-R10011B-005N	Silica	4 x 8 cm	150 μm	None	50	
TLP-R10011B-117	Silica	2.5 x 7.5 cm	200 μm	F254	200	
TLP-R10011B-323	Silica	20 x 20 cm	200 μm	F254	25	
TLP-R10011B-323N	Silica	20 x 20 cm	200 μm	None	25	
TLP-R10014B-0115	Silica	5 x 6.7 cm	200 μm	F254	50	
TLP-R10014BK-0116	Silica, optimized for KMnO4 revelation	3.3 x 6.6 cm	200 μm	F254	50	

#### **TLC Accessories**

4	TLC Accessories		
PN	Accessory		Qty / Box
AUT-0161	Rectangular TLC Developing Chamber		1
AUT-0161B	Replacement Lid for Rectangular Developing Chamber		1
AUT-0162	TLC Adsorbent Scraper		1
AUT-0163	TLC Spotting Capillaries		300
AUT-0182	Cutter for glass-baked TLC Plates (up to 20 x 20 cm)		1
AUT-0183	Replacement Scriber for TLC Plate Cutter		1
AUT-1182	TLC Plate Pencil Glass Cutter	E CO LOCAL	1

40

# **Thin Layer Chromatography Practical Guide**

# **Select a Stationary Phase**

As almost 80 % of all separations can be performed using silica gel plates, it is suggested to try using this coating first. However, for acid sensitive compounds, alumina is probably a better choice (useful for amine purification). If you are working with highly polar compounds, reversed-phase mode is more suitable.

### Select a Mobile Phase (Solvent Systems)

The selection of the mobile phase (also called solvent system or eluent) is perhaps the most important parameter to achieve efficient thin-layer chromatography separation. It is based on the compound's solubility with the solvent and the difference in the affinity for the mobile phase versus the stationary adsorbent (silica or alumina).

In normal phase chromatography, where non-polar solvents such as hexane or pentane are used, non-polar compounds will move up the plate while most polar compounds will stay on the baseline. Inversely, polar solvents will allow polar compounds to move off the origin. The most suitable solvent system is the one that moves all components off the baseline with Rf values between 0.15 and 0.85 (ideally, close to 0.2 - 0.4).

For most applications, a common solvent system to start with is EtOAc / Hexane (1:1). Varying the ratio can have a pronounced effect on the Rf. If it is not working, then try: MeOH / DCM (2:8 - 10:90); or toluene with acetone, EtOAc, or DCM.

Remember: in normal phases, to increase the compound's Rf, increase the polarity of the mobile phase; increase the ratio of the polar solvent or choose another solvent. Inversely, to decrease Rf, decrease the polarity of the eluent.

#### Rules of thumb

- Standard compounds (most popular solvent system): 10 50 % EtOAc / Hexane
- Polar compounds: 100 % EtOAc or 5 10 % MeOH / DCM
- Non-polar compounds: 5 % EtOAc (or ether) / Hexane or 100 % Hexane
- For basic compounds: (amine or nitrogen containing), it could be useful or required to add a small quantity of triethylamine ( $Et_2N$ ) to the solvent mixture (0.1 - 2.0 % but typical quantity is 0.1 %) or 1 - 10 % ammonia (NH<sub>2</sub>) in MeOH / DCM.
- For acidic compounds: it could be useful to add acetic (AcOH) or formic acid (FA) to the solvent mixture (0.1 2.0 %).

#### Reversed-phase mode

Solutions for Chromatography and Purification

In reversed-phase chromatography, the typical solvent systems are:

- Mixtures of water or aqueous buffers and water miscible organic solvents such as acetonitrile (ACN), methanol and tetrahydrofuran (THF). Other solvents can be used such as ethanol (EtOH) and isopropanol (IPA).
- MeOH, to improve peak shape in flash chromatography, 0.1 % of acetic, formic or trifluoroacetic acid (TFA) can be added to the solvent system.





# **TLC Preparation & Interpretation**

### **TLC Plate Preparation**

Using a pencil, lightly draw a straight-line parallel to the width of the plate at about 1 cm from the base end of the plate. Sample application will be done on this line called baseline or origin.

Note: never use a pen because ink can move with some solvents used as eluent.

### **Sample Preparation**

Thorough sample preparation is a prerequisite for an optimal and efficient TLC separation. Typical sample preparation processes could consist in a sample crushing, filtration, extraction or concentration of the product of interest.


# **Sample Application**

Sample preparation will differ depending on the nature of the plate (*analytical or preparative*). For analytical plates, because thin layer chromatography is extremely sensitive, it is really important to apply a small quantity using a glass capillary (*or a micro pipette*) to get optimal resolution. For preparative plates, apply a series of small adjacent spots to form a band or a streak using a glass capillary (*or a microliter syringe*). In both cases, a spotting guide can be used to facilitate sample application.

# Co-spotting

For analytical chromatography, co-spotting is frequently used for similar polarity products.

This consists to apply on the same spot the starting material and reaction mixture, as shown by the image below.





### **TLC Plate Development**

The most commonly used method to perform thin layer chromatography separation is to place vertically the TLC plate inside a sealed developing chamber to ensure solvent saturation. Place approximately 0.5 cm of the suitable solvent system inside the chamber. Slowly place the TLC inside the chamber and allow the eluent to travel up the plate until it gets to 1 cm from the top of the plate. Immediately remove the plate and draw a line along the solvent front.

**Note**: for optimal solvent saturation, a filter paper can be added inside the TLC chamber. This also prevents eluent evaporation. The solvent level needs to be below the baseline; otherwise the spots will be dissolved.

### **TLC Plate Visualization**

If components of the reaction are colored, no visualization method is required (*spots can be seen directly on the silica layer*). However, most of the time it is not the case, therefore one of the methods described below should be used to reveal the spots.

#### Non-destructive methods

As a general visualization procedure, before treating the TLC plate with any destructive methods, UV-active compounds can be viewed under an ultraviolet lamp (*usually for polyconjugated compounds like benzophenones and anthracenes*). Furthermore, an iodine chamber can be useful for thiols, phosphines and alkenes but it works in about 50 % of cases for alkanes. It is recommended to circle the spots with a pencil on the TLC plate prior to visualization by destructive methods.

#### **Destructive methods**

For compounds that are not UV-active, there are several varieties of stains that can be used depending on the nature of the compound of interest. To use a stain, simply dip the TLC plate into the staining solution as quickly as possible, and then immediately absorb the excess stain with paper and heat carefully with a heat gun or on a hot plate at 110°C until spots are revealed.

# **Chromatogram Interpretation**

### Retention factor (Rf) definition

Retention factor analysis is used to evaluate if the solvent system is adequate. Rf is defined as the distance traveled by the compound divided by the distance traveled by the solvent front. This means: the larger the Rf value of a compound, the larger is the distance traveled by the compound. In other words, when comparing Rf values of various compounds under identical chromatography conditions, the compound with the larger Rf is less polar because it interacts less strongly with the polar adsorbent on the plate.

**Remember**, a good solvent system is one that moves all components off the baseline with Rf values between 0.15 and 0.85 (*ideal Rf is 0.2 - 0.4*). Otherwise, when possible, it is preferable to chose another solvent system.

Retention factor (Rf) =  $\frac{\text{distance traveled by the compound}}{\text{distance traveled by the solvent front}}$ 

Rf calculation based on the example shown here: Rf = 4.0 cm / 5.5 cm = 0.73


#### **Prediction of Column Volumes (CV)**

TLC data can be used to predict column elution based on the relationship between the retention factor and the column volume. CV is the number of column volumes required to elute the component from the column regardless of column dimensions [(bed volume) - (volume of packing)].

CV = 1 / Rf &  $\Delta CV = 1 / Rf_1 - 1 / Rf_2$ 

The greater the  $\Delta$ CV, the greater will be the separation and resolution between the spots (easier separation).

A bigger  $\Delta$ CV will therefore allow more sample to be loaded onto the column.







# **TLC Visualization Methods**

Described below are the most frequently used TLC visualization methods (also called stains) in alphabetical order.

	Stains for T	Thin Layer Chromatography	
Name	Visualization of	Stain Recipe	Comments
p-Anisaldehyde #1	Universal stain	Prepare stain as follows	Visualization Colors
	Good for nucleophiles and	2 mL of glacial acetic acid	Spots: Various colors
	oxygenated compounds	5 mL of p-anisaldehyde	BG: Orange to pink
		7 mL of conc. sulfuric acid	Appropriate Storage
		• 185 mL of 95 % ethanol	Aluminum wrapped at 0 <sup>o</sup>
		Tip: add dropwise the acid at the end and stir vigorously.	
<b>Note</b> : Tends to be insensitive to a	     Ikenes, alkynes and aromatic compo	 unds unless other functional groups are present.	
p-Anisaldehyde #2	Acronycine	Prepare stain as follows	Visualization Colors
	Cineoles	1 mL p-anisaldehyde	Spots: Various colors
	Terpenes	10 mL perchloric acid	BG: Orange to pink
		20 mL acetone	Appropriate Storage
		80 mL water	Aluminum wrapped at 0°
Bromocresol Green	Acidic groups (pK <sub>a</sub> < 5)	Prepare stain as follows	Visualization Colors
	Carboxylic acids	0.04 g of bromocresol green	Spots: Yellow to green
		• 100 mL of 95 % ethanol	BG: Blue
		0.1 M solution of sodium hydroxide	Appropriate Storage
		Tip: add the base slowly at the end until the solution turns	Aluminum wrapped at 0
		pale blue.	Heating NOT required
Cerium Molybdate	Universal stain	Prepare stain as follows	Visualization Colors
(CAM or Hanessian's Stain)	Good for peptides	12 g of ammonium molybdate	Spots: Blue
		0.5 g of ceric ammonium molybdate	BG: White
		15 mL of conc. sulfuric acid	Appropriate Storage
		235 mL of water	Aluminum wrapped
Note: Highly sensitive stain; very	low concentration of product may app	pear as a significant impurity.	
Cerium Sulfate	Difficultly stainable compounds	Prepare stain as follows	Visualization Colors
$(Ce(SO_4)_2)$		15 % aqueous sulfuric acid saturated with ceric sulfate	Spots: Black
			BG: Yellow to white
Chromic Acid	Difficultly stainable compounds	Prepare stain as follows	
		2.5 g of potassium chromate	
		• 100 mL of 20 % sulfuric acid in water	
Cobalt Chloride	Universal stain	Prepare stain as follows	Visualization Colors
(CoCl <sub>2</sub> )	Used in conjunction with PMA	2 g of cobalt chloride	Spots: Various colors
-	when this one is not effective	• 100 mL of water	BG: Pink
	enough	10 mL of conc. sulfuric acid	Heating NOT required
		<b>Tip</b> : simply dip PMA treated plate in CoCl <sub>2</sub> solution.	
p-Dimethylamino-benzaldehyde	Amines	Prepare stain as follows	Visualization Colors
	Ladatas	0.5 g of p-dimethylamino-benzaldehyde	Spots: Blue
(PDAB or Ehrlich's Reagent)	Indoles	0.0 g of p difficulty armino beneauderly de	-
(PDAB or Ehrlich's Reagent)	Indoles	10 mL of conc. hydrochloric acid	BG: White

Abbreviation: BG stands for "background".



$\sim$	Stains	for Thin Layer Chromatography	
Name	Visualization of	Stain Recipe	Comments
2,4-Dinitrophenyl-hydrazine	Aldehydes	Prepare stain as follows	Visualization Colors
(DNP)	Ketones	12 g of 2,4-dinitrophenylhydrazine	Spots: Yellow to red
		60 mL of conc. sulfuric acid	BG: Light orange
		80 mL of water	DO NOT HEAT dipped plate
		200 mL of 95 % ethanol	
Dragendorff Reagent	Nitrogenous compounds	Prepare stain as follows:	Visualization Colors
	(alkaloids, amines, organics bases, etc.)	Solution A	Spots: Orange to red
	Phenols	1.7 g of bismuth nitrate	BG: Yellow
	Prieriois	80 mL of water	Appropriate Storage
		20 mL of acetic acid	Aluminum wrapped
		Solution B	Stain Shelf-Life
		40 g of potassium iodide	One or two weeks
		100 mL of water	Solutions A and B are
		<b>Tip:</b> mix 5 mL of each solution A and B to a solution of	long term storable
		20 mL of acetic acid in 70 mL of water.	DO NOT HEAT dipped plate
Ferric Chloride	Phenols	Prepare stain as follows	Visualization Colors
(FeCl <sub>2</sub> )		2 g of ferric chloride	Spots: Red
. 3/		102 mL of 0.5 N hydrochloric acid	BG: Yellow
ladina	Unacturated and aromatic	Dranara atain aa fallawa	Vigualization Colors
lodine	Unsaturated and aromatic compounds	Prepare stain as follows	Visualization Colors
		lodine crystals in an amber bottle	Spots: Dark brown     BG: Light brown
	1	I	1 • B(¬: 1 10111 D10W1)
	I		Bo. Light Slown
Note: lodine stain can be remov	l ved by heating.		De. Light Stewn
	ved by heating.  Universal stain	Prepare stain as follows	Visualization Colors
Note: Iodine stain can be remov  Morin Hydrate (Hydroxy Flavone)	T	Prepare stain as follows  • 0.1 % of morin hydrate in methanol	· · · · · · · · · · · · · · · · · · ·
Morin Hydrate	Universal stain		Visualization Colors
Morin Hydrate	Universal stain	0.1 % of morin hydrate in methanol	Visualization Colors  • Spots: Various colors
Morin Hydrate (Hydroxy Flavone) Ninhydrin	Universal stain Fluorescently active  Amino acids	0.1 % of morin hydrate in methanol  Prepare stain as follows	Visualization Colors  • Spots: Various colors  • BG: White  Visualization Colors
Morin Hydrate (Hydroxy Flavone)	Universal stain Fluorescently active	0.1 % of morin hydrate in methanol	Visualization Colors  • Spots: Various colors  • BG: White
Morin Hydrate (Hydroxy Flavone) Ninhydrin	Universal stain Fluorescently active  Amino acids Amino sugars	• 0.1 % of morin hydrate in methanol  Prepare stain as follows • 1.5 g of ninhydrin	Visualization Colors  • Spots: Various colors  • BG: White  Visualization Colors  • Spots: Various colors
Morin Hydrate (Hydroxy Flavone) Ninhydrin (Indanetrione Hydrate)	Universal stain Fluorescently active  Amino acids Amino sugars Amines	Prepare stain as follows  1.5 g of ninhydrin  3 mL acetic acid  100 mL of n-butanol	Visualization Colors  • Spots: Various colors  • BG: White  Visualization Colors  • Spots: Various colors  • BG: White
Morin Hydrate (Hydroxy Flavone)  Ninhydrin (Indanetrione Hydrate)  Phosphomolybdic Acid	Universal stain Fluorescently active  Amino acids Amino sugars Amines  Universal stain	Prepare stain as follows  1.5 g of ninhydrin  3 mL acetic acid  100 mL of n-butanol  Prepare stain as follows	Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Various colors  BG: White  Visualization Colors
Morin Hydrate (Hydroxy Flavone)  Ninhydrin (Indanetrione Hydrate)  Phosphomolybdic Acid	Universal stain Fluorescently active  Amino acids Amino sugars Amines	Prepare stain as follows  1.5 g of ninhydrin  3 mL acetic acid  100 mL of n-butanol  Prepare stain as follows  10 % of PMA solution in ethanol	Visualization Colors  • Spots: Various colors  • BG: White  Visualization Colors  • Spots: Various colors  • BG: White  Visualization Colors  • Spots: Dark green to black
Morin Hydrate (Hydroxy Flavone)  Ninhydrin (Indanetrione Hydrate)  Phosphomolybdic Acid (PMA)	Universal stain Fluorescently active  Amino acids Amino sugars Amines  Universal stain Very effective against diluted sample	Prepare stain as follows  1.5 g of ninhydrin  3 mL acetic acid  100 mL of n-butanol  Prepare stain as follows  10 % of PMA solution in ethanol  or 10 g of PMA in 100 mL of ethanol	Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Dark green to black  BG: Light green
Morin Hydrate (Hydroxy Flavone)  Ninhydrin (Indanetrione Hydrate)  Phosphomolybdic Acid (PMA)  Potassium Permanganate*	Universal stain Fluorescently active  Amino acids Amino sugars Amines  Universal stain Very effective against diluted sample  Olefins	Prepare stain as follows  1.5 g of ninhydrin  3 mL acetic acid  100 mL of n-butanol  Prepare stain as follows  10 % of PMA solution in ethanol  or 10 g of PMA in 100 mL of ethanol  Prepare stain as follows  Prepare stain as follows	Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  BG: White  Visualization Colors  Spots: Dark green to black  BG: Light green  Visualization Colors
Morin Hydrate (Hydroxy Flavone)  Ninhydrin (Indanetrione Hydrate)  Phosphomolybdic Acid (PMA)  Potassium Permanganate*	Universal stain Fluorescently active  Amino acids Amino sugars Amines  Universal stain Very effective against diluted sample  Olefins Readily oxidized groups	Prepare stain as follows  1.5 g of ninhydrin  3 mL acetic acid  100 mL of n-butanol  Prepare stain as follows  10 % of PMA solution in ethanol  or 10 g of PMA in 100 mL of ethanol  Prepare stain as follows  1.5 g of potassium permanganate	Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Dark green to black  BG: Light green  Visualization Colors  Spots: Yellow to light brown
Morin Hydrate (Hydroxy Flavone)  Ninhydrin (Indanetrione Hydrate)  Phosphomolybdic Acid (PMA)	Universal stain Fluorescently active  Amino acids Amino sugars Amines  Universal stain Very effective against diluted sample  Olefins Readily oxidized groups Alcohols, aldehydes,	Prepare stain as follows  1.5 g of ninhydrin  3 mL acetic acid  100 mL of n-butanol  Prepare stain as follows  10 % of PMA solution in ethanol  or 10 g of PMA in 100 mL of ethanol  Prepare stain as follows  1.5 g of potassium permanganate  10 g of potassium carbonate	Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Dark green to black  BG: Light green  Visualization Colors  Spots: Yellow to light brown  BG: Purple to pink
Morin Hydrate (Hydroxy Flavone)  Ninhydrin (Indanetrione Hydrate)  Phosphomolybdic Acid (PMA)  Potassium Permanganate*	Universal stain Fluorescently active  Amino acids Amino sugars Amines  Universal stain Very effective against diluted sample  Olefins Readily oxidized groups	Prepare stain as follows  1.5 g of ninhydrin  3 mL acetic acid  100 mL of n-butanol  Prepare stain as follows  10 % of PMA solution in ethanol  or 10 g of PMA in 100 mL of ethanol  Prepare stain as follows  1.5 g of potassium permanganate  10 g of potassium carbonate  1.25 mL of 10 % sodium hydroxide	Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Dark green to black  BG: Light green  Visualization Colors  Spots: Yellow to light brown
Morin Hydrate (Hydroxy Flavone)  Ninhydrin (Indanetrione Hydrate)  Phosphomolybdic Acid (PMA)  Potassium Permanganate*	Universal stain Fluorescently active  Amino acids Amino sugars Amines  Universal stain Very effective against diluted sample  Olefins Readily oxidized groups Alcohols, aldehydes,	Prepare stain as follows  1.5 g of ninhydrin  3 mL acetic acid  100 mL of n-butanol  Prepare stain as follows  10 % of PMA solution in ethanol  or 10 g of PMA in 100 mL of ethanol  Prepare stain as follows  1.5 g of potassium permanganate  10 g of potassium carbonate	Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Dark green to black  BG: Light green  Visualization Colors  Spots: Yellow to light brown  BG: Purple to pink
Morin Hydrate (Hydroxy Flavone)  Ninhydrin (Indanetrione Hydrate)  Phosphomolybdic Acid (PMA)  Potassium Permanganate* (KMnO <sub>4</sub> )	Universal stain Fluorescently active  Amino acids Amino sugars Amines  Universal stain Very effective against diluted sample  Olefins Readily oxidized groups Alcohols, aldehydes, alkenes, alkynes, etc.	Prepare stain as follows  1.5 g of ninhydrin  3 mL acetic acid  100 mL of n-butanol  Prepare stain as follows  10 % of PMA solution in ethanol  or 10 g of PMA in 100 mL of ethanol  Prepare stain as follows  1.5 g of potassium permanganate  10 g of potassium carbonate  1.25 mL of 10 % sodium hydroxide	Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Dark green to black  BG: Light green  Visualization Colors  Spots: Yellow to light brown  BG: Purple to pink  Stain Shelf-Life
Morin Hydrate (Hydroxy Flavone)  Ninhydrin (Indanetrione Hydrate)  Phosphomolybdic Acid (PMA)  Potassium Permanganate* (KMnO <sub>4</sub> )	Universal stain Fluorescently active  Amino acids Amino sugars Amines  Universal stain Very effective against diluted sample  Olefins Readily oxidized groups Alcohols, aldehydes, alkenes, alkynes, etc.	Prepare stain as follows  1.5 g of ninhydrin  3 mL acetic acid  100 mL of n-butanol  Prepare stain as follows  10 % of PMA solution in ethanol  or 10 g of PMA in 100 mL of ethanol  Prepare stain as follows  1.5 g of potassium permanganate  10 g of potassium carbonate  1.25 mL of 10 % sodium hydroxide  200 mL of water	Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Dark green to black  BG: Light green  Visualization Colors  Spots: Yellow to light brown  BG: Purple to pink  Stain Shelf-Life
Morin Hydrate (Hydroxy Flavone)  Ninhydrin (Indanetrione Hydrate)  Phosphomolybdic Acid (PMA)  Potassium Permanganate* (KMnO <sub>4</sub> )	Universal stain Fluorescently active  Amino acids Amino sugars Amines  Universal stain Very effective against diluted sample  Olefins Readily oxidized groups Alcohols, aldehydes, alkenes, alkynes, etc.	Prepare stain as follows  1.5 g of ninhydrin  3 mL acetic acid  100 mL of n-butanol  Prepare stain as follows  10 % of PMA solution in ethanol  or 10 g of PMA in 100 mL of ethanol  Prepare stain as follows  1.5 g of potassium permanganate  10 g of potassium carbonate  1.25 mL of 10 % sodium hydroxide  200 mL of water  mercaptans groups when gently heated.	Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Dark green to black  BG: Light green  Visualization Colors  Spots: Yellow to light brown  BG: Purple to pink  Stain Shelf-Life  Three months
Morin Hydrate (Hydroxy Flavone)  Ninhydrin (Indanetrione Hydrate)  Phosphomolybdic Acid (PMA)  Potassium Permanganate* (KMnO <sub>4</sub> )	Universal stain Fluorescently active  Amino acids Amino sugars Amines  Universal stain Very effective against diluted sample  Olefins Readily oxidized groups Alcohols, aldehydes, alkenes, alkynes, etc.  of alcohols, amines, sulfides and municipal stain	Prepare stain as follows  1.5 g of ninhydrin  3 mL acetic acid  100 mL of n-butanol  Prepare stain as follows  10 % of PMA solution in ethanol  or 10 g of PMA in 100 mL of ethanol  Prepare stain as follows  1.5 g of potassium permanganate  10 g of potassium carbonate  1.25 mL of 10 % sodium hydroxide  200 mL of water  prepare stain as follows  Prepare stain as follows  Prepare stain as follows	Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Various colors  BG: White  Visualization Colors  Spots: Dark green to black  BG: Light green  Visualization Colors  Spots: Yellow to light brown  BG: Purple to pink  Stain Shelf-Life  Three months

Otherwise mentioned, stains are long-term stable when stored in a tightly-closed container to prevent solvent evaporation



# Silia Plate TLC Troubleshooting

### Streaking or elongated spot rather than a defined spot?

#### **Possible Solutions**

- Sample was overloaded: run the TLC again using a more diluted solution of your sample.
- In presence of a base sensitive compound: try to add acetic or formic acid to the eluent (0.1 2.0 %).
- In presence of an acid sensitive compound: try to add triethylamine to the eluent (0.1 2.0 %) or 1 10 % ammonia in MeOH / DCM. If it is not working use Alumina as TLC coating.
- In presence of too highly polar compounds: try using a specialized silica TLC plate like reversed-phase (C18 for example).

### Unable to see any spots on the TLC?

#### Possible Solutions

- If you have not been able to visualize any spots on your TLC using UV light, try another method; maybe your compound is not UV-
- Maybe your sample is too diluted. Try to apply several times your sample on the same spot (do not forget to dry solvent between each application for optimal results) or to concentrate your solution.
- · Make sure the solvent level inside the tank is lower than the spotting line to avoid sample dissolution by the eluent.

# How to monitor a reaction in presence of similar Rfs for both starting materials and product of interest?

#### **Possible Solutions**

- · Try the co-spotting method.
- Try to visualize the plate using anisaldehyde or molybdene. Spot color or brightness differ for two compounds when using these stains.
- If none of the two previous solutions work, change solvent systems (use another class of solvent).

Tips: in chromatography, there are three classes of solvent systems providing significantly different results:

- 1. Mixture of polar / hydrocarbon solvents (i.e.: EtOAc / Hexane; Ether / Petroleum ether).
- 2. Mixture of polar / dichloromethane solvents (examples of polar solvent: Ether, EtOAc, MeOH).
- 3. Mixture of polar / benzene (or toluene) solvents (examples of polar solvent: Ether, EtOAc, MeOH).

# Compounds stay too close to the baseline or solvent front.

#### **Possible Solutions**

- · Too close to the baseline: your eluent is not polar enough; increase the proportion of polar solvent in the same solvent system or chose a more polar solvent.
- Too close to the solvent front: inversely, your eluent is too polar; decrease the proportion of polar solvent in the same solvent system or chose a less polar solvent.

### **Resource Center**

#### **Download our Posters**







Choosing your metal & organic scavenging solution

SiliCycle's most popular functionalized silicas

SiliCycle packing formats in a nutshell

•

Functionalized silicas and reference information

Take a Look at some of our Multimedia Contents



Introduction to metal and organic Scavengers



Metal scavenging using bulk SiliaMetS functionalized silica



How to calculate the amount of What are the parameters that influence scavenging efficiency? scavenger needed

Get a Copy of our E-Books



E-PAK flow purification cartridges



Scale-up impurity scavenging with E-PAK



E-PAK cartridge housings, from lab to commercial scale



See how easy it is working with E-PAK



Flash separation of dye mixture with SiliaSep Premium



How does flash chromatography work?



**Understanding Column** Volume



between retention factor and column volume



The 5 steps of a solid phase extraction (SPE)



Understanding particle size distribution - D50, D90 and D10



A collection of various case studies and application notes using scavengers



A collection of scale-up case studies and application notes using scavengers



What pH range is suitable for functionalized silica?



capacity of preparative TLC plates?







#### We Redefine Purity

High Value Silica-based and Specialty Products for Organic Chemistry, Chromatography, and Purification

#### **DISCOVER AND DOWNLOAD OUR BROCHURES**

#### METAL AND ORGANIC SCAVENGING

SiliaMetS® - Metal Scavengers

SiliaBond® - Organic Scavengers

E-PAK® - Fixed Bed Flow-Through Purification Cartridges



#### **CHROMATOGRAPHY AND PURIFICATION**

SiliaFlash® - Irregular Silica Gels | SiliaSphere™ PC - Spherical Silica Gels

SiliaBond® - Chromatographic Phases


SiliaSep™ - Flash Cartridges | SiliaPlate™ - TLC Plates



#### SAMPLE PREPARATION

SiliaPrep™ - Silica-based SPE Cartridges and Well Plates


SiliaPrepX™ - Polymeric SPE Cartridges and Well Plates



#### ANALYTICAL AND PREPARATIVE CHROMATOGRAPHY

SiliaSphere™ - Spherical Silica Gels

SiliaChrom® - HPLC Columns



#### **ORGANIC SYNTHESIS**

SiliaBond® - Reagents and Oxidants

SiliaCat® - Heterogeneous Catalysts



#### **PEPTIDE SYNTHESIS**

Peptide Synthesis and Purification Solutions Amine Free Basing and TFA Removal



#### **R&D SERVICES**

Metal and Organic Scavenging Screenings | Organic Synthesis Chromatography and Purification | Material Science Method Development, Optimization, and Transfer



#### **Technical Support**

At SiliCycle, we are committed to providing the best technical support possible.

Our worldwide Technical Support Group of highly qualified M. Sc., Ph. D. Chemists and Engineers will answer your questions and provide solutions to your most advanced chemistry and purification needs. Contact us at support@silicycle.com or call us.



SiliCycle Inc.

2500, Parc-Technologique Blvd, Quebec City (Quebec) G1P 4S6 CANADA 🎞 🛂

europe@silicycle.com

Phone: +1 418.874.0054

Fax: **+1 418.874.0355** 

Toll Free **+1 877.745.4292** (North America only)

Email: info@silicycle.com




Overseas Offices

SiliCycle Europe

+33 1 87 65 06 57





Follow us: in 🔰 🧗 You Tube

