



A Practical Cleavage of Acetals Using SiliaMetS[®] Thiol

Removing protecting groups from a synthesis pathway can prove difficult. In order to overcome this challenge, a method using SiliaMetS[®] Thiol was developed for the deprotection of various aromatic and aliphatic acetals to their corresponding catechol or diol derivatives.

This interesting process has many advantages, including an overall simplified procedure compared to the traditional deprotection methods. To start it is a mild procedure, is amenable to scale-up, does not require inert atmosphere and can be carried out at ambient temperature. Both clean conversions and good chemoselectivity can be observed. The thiol stench was avoided, and it allows the use of solvent independent thermally stable reagents. In addition, this process allows the removal of toxic or difficult reagents and by-products by a simple filtration. Finally, in many cases, purification by chromatography was unnecessary which reduces the amount of solid waste while not impacting yields adversely.




Table 1: Yields of Catechol 2 in different solvents and acids

Yields of Catechol 2 in different solvents and acids				
Entry	Si-Thiol (equiv)	Lewis acid (equiv)	Solvent	Yield (%)
1	1.0	BF ₃ · Et ₂ O (2.0)	Dichloromethane	75
2	1.0	BF ₃ · Et ₂ O (1.0)	Dichloromethane	88
3	5.0	BF ₃ · Et ₂ O (1.0)	Dichloromethane	93
4	2.0	BF ₃ · Et ₂ O (1.0)	Dichloromethane	>99
5	2.0	FeCl ₃ (1.0)	Dichloromethane	99
6	2.0	AlCl ₃ (1.0)	Dichloromethane	99
7	2.0	ZnCl ₂ (1.0)	Dichloromethane	99
8	2.0	p-TsOH (1.0)	Dichloromethane	88
9	2.0	BF ₃ · Et ₂ O (1.0)	Acetonitrile	>99
10	2.0	BF ₃ · Et ₂ O (1.0)	Dichloromethane	90
11	2.0	BF ₃ · Et ₂ O (1.0)	tert-Butyl methyl ether	90
12	2.0	BF ₃ · Et ₂ O (1.0)	Diethyl ether	74
13	2.0	BF ₃ · Et ₂ O (1.0)	Acetone	25
14	2.0	BF ₃ · Et ₂ O (1.0)	N,N-Dimethyl-formamide	Mixture of products
15	2.0	BF ₃ · Et ₂ O (1.0)	N,N-Dimethyl-acetamide	Mixture of products
16	2.0	BF ₃ · Et ₂ O (1.0)	Methanol	No reaction

This method was also investigated by examination of an array of structurally diverse acetals and it gave promising results with good chemoselectivity being displayed.

De Léséleuc, M. et al. Eur. J. of Org. Chem., 2019, 44, 7389 - 7393
Prometic Biosciences Inc.

LEARN MORE

about SiliaMetS[®] Thiol Metal Scavenger (R51030B) on our website:
<https://www.silicycle.com/products/metal-scavengers/portfolio/r51030b>