

# New Generation of HPLC Silica-Based C18 Column for Both Highly Acidic and Basic pH:

## The SiliaChrom® SB C18 and SiliaChrom® XT C18 Fidelity

Charles Lévesque, Vincent Bédard and Francois Béland

SiliCycle® Inc. (www.silicycle.com) 2500 Parc-Technologique Blvd, Quebec City, Quebec, G1P 4S6, Canada

\*(charleslevesque@silicycle.com)

Si

### Introduction

Acidic mobile phases have found widespread applications in the reversed-phase HPLC separation of many important pharmaceutical and environmental compounds. Analytes such as pharmaceuticals and biomolecules often show peak shape, retention and selectivity changes when the mobile phase pH is changed from neutral to acidic pH (pH=1.0). In fact, lowering the pH helps to suppress silanol interactions between basic solutes and the residual surface silanols, thus resulting in less tailing and better retention for acidic compounds (pKa lower than 2), thus increasing retention and altering selectivity (figure 1). For example, hydrophobic peptides and protein separations use TFA for solubility. Additionally, basic analytes are preferably positively charged in the HPLC eluent for positive-mode electrospray LC-MS applications thus requiring low pH mobile phases and high column robustness.

Also is well known in reversed phase chromatography that the uncharged form of a compound is best suited for troubleless HPLC (figure 1). Most basic drugs have pK<sub>a</sub> values of around 9.5. To keep basic substances uncharged they must be kept in an environment with pH higher than its pK<sub>a</sub> value of the analyte. Today, HPLC columns with silica based stationary phase for acidic and basic mobile phases are important solution for end-users in HPLC.




Figure 1: Effect of pH on the retention acids and bases in reversed phase HPLC

### HPLC Column Evaluation SiliaChrom SB C18 and SiliaChrom XT C18 Fidelity

#### SiliaChrom SB C18

The surface of SiliaChrom SB C18 is treated with an organic form of silicon to increase the number of silanol groups on the surface. After this, the surface is bonded with two types of octadecyl silanes. One of them has a protecting group that shields the area under it and protects the surface from an acid attack from the mobile phase (figure 2).




Figure 2 SiliaChrom SB C18 phase

#### SiliaChrom XT C18 Fidelity

The key is to have a hybrid particle or hybrid surface to reduce the solubility of silica at high pH. In fact, the SiliaChrom XT is silica coated with a prepolymer of ethyltriethoxysilane/tetraethoxysilane, followed by special thermo treatment to get a rigid surface that is less soluble than silica itself at high pH (figure 3)




Figure 3 SiliaChrom XT C18 Fidelity phase

### Stability Study

The new SiliaChrom SB and XT C18 5µm phase were packed in-house. The column configuration was 4.6mm x 150mm. The low and high pH stability test were performed using a Thermo Surveyor HPLC with variable UV detector monitoring absorbance at 270 nm, using a flow rate of 1.00 mL/min with the column heated at 35°C. A sample containing uracil (void marker), phenol, toluene and ethylbenzene was injected to get K' between 1.0 and 1.5 in order to be very sensitive to minute variations of the stationary phase. We have monitored all injections every 6 hours for a total of 120 hours. The flow was kept constant during 120 hours (1.00 mL/min, 35°C) with appropriate aggressive mobile phase:

#### Low pH conditions

ACN:Water:TFA (59:39:2) pH=1.00

UV 270 nm

T=23.0°C

1.00 mL/min

SiliaChrom SB C18, 5µm, 4.6x150mm (SiliCycle P/N: HPC-1105-15N)

#### High pH conditions

ACN:Water:TEA (54.9:44.9:0.2) pH=11.50


UV 270 nm

T=23.0°C

1.00 mL/min


SiliaChrom XT C18 Fidelity, 5µm, 4.6x150mm (SiliCycle P/N: HPC-11005-15N)

#### Test Mix:



#### Results

The chromatograms illustrate separations where the retention time for each analyte is constant. The Figure 2 shows no modification of retention time from the beginning (0 hour) and 120 hours later with the flow kept constant with the low or high pH mobile phase



#### Conclusion

This poster demonstrates clearly that the SiliaChrom SB C18 and SiliaChrom XT C18 Fidelity are very stable for aggressive mobile phases (low and high pH). The beauty is the fact that material for each one still to be silica based in order to keep the powerful efficiency in HPLC. We decided to stop the acquisition after 120 hours but the stability limit of those HPLC columns is not well known. We presume no degradation at all, for at least many months as typical HPLC column using buffer with pH range 2-8. Also we did not observed any tailing effect for toluene as XTerra® shown. It is well known that small molecules have some tailing with hybrid packing (hole trap). SiliaChrom SB C18 and XT C18 Fidelity are the best HPLC columns for method development especially using pH as a tool to have best HPLC analysis.